DISCLAIMER FOR COMMUNICATION INTERFACES

Fluke provides communications interfaces ("as provided in the relevant communications interface document") for development and use by its customers ("Users") for their own use and within User's controlled environment. Fluke uses and has verified the functionality of these communications interfaces in accordance with its own in-house performance standards using accepted test procedures.

Except for liability which cannot be excluded by law, Fluke shall have no responsibility for User's development or use of these communications interfaces nor for any loss, corruption or limitation of access to the communications interfaces. Fluke is not responsible for any trouble shooting nor is it responsible for any damages incurred to any device used for testing (whether a Fluke device or third-party device), the equipment being tested or any person. Fluke will not be responsible for any injuries sustained due to unauthorized equipment modifications.

These communications interfaces are provided to User as-is and provided without warranty of any kind, whether statutory, written, oral, express or implied (including any warranties of merchantability or fitness for a particular purpose or any warranties arising from course of dealing or usage of trade). Fluke does not warrant that the communications interfaces will be delivered free of any interruptions, delays, omissions or errors or in a secure manner or that any faults or trouble shooting will be corrected.

Follow any other instructions in the relevant communications interface document and do not touch the Fluke device while deploying the communications interfaces as this may result in electrical shock hazards or improper operation.

No responsibility is assumed by Fluke for the use or reliability of communications interfaces that are not supplied by Fluke.

FLUKE Biomedical Impulse 6000D/7000DP Communications Interface

Version 2.4 4/28/2016

INTRODUCTION

The Impulse 6000D and Impulse 7000DP can be controlled remotely through its Computer Port: a USB Interface port that looks like a COM port to a personal computer (PC).

You can control the Impulse by sending commands to it and receiving responses, including test data, through the COM port.

OPERATING SYSTEM REQUIREMENT

Fluke supports connecting the Impulse to a PC running Windows XP, Vista, or a later version.

VIRTUAL COM PORT

The Impulse USB port is built from an integrated circuit (IC) device that is commonly used inside adapter cables that convert USB to RS232. When this device is connected to a PC it looks like a COM port to the PC. When Windows enumerates the device it assigns a COM port number to it. It is called a virtual COM port (VCP).

The IC is an FT232R from the FTDI company. It is compatible with the USB Version 2.0 Full Speed specification.

The USB port may reside inside the Impulse, but the PC acts like it now has an additional COM port and that COM port is connected to an RS232 serially controlled instrument.

Once your Impulse is connected and the COM port is enumerated, you can control it by sending remote commands to the COM port and receiving responses.

COMPUTER CABLE CONNECTION

The Impulse Computer Port is a USB Device Port (peripheral) with a Type B square connector. It connects to a PC USB Controller Port that has a Type A rectangular connector.

Connect the Impulse to your PC with a standard USB Type A to Type B cable such as the one supplied with the Impulse.

WINDOWS SOFTWARE DRIVER

Versions of Windows XP, Vista, and later, include a software driver for FTDI USB Serial Converters, including the FT232R. The USB ID numbers are: VID 0403 and PID 6001.

When you connect the Impulse to your PC for the first time, Windows should recognize and register your Impulse as a USB Serial Converter and USB Serial Port (COMx).

DEVICE MANAGER

Run Device Manager to check the status of the Impulse COM port. When viewing by Type, your Impulse shows up in two places:

- Universal Serial Bus controllers / USB Serial Converter.
- Ports (COM & LPT) / USB Serial Port (COMx).

If you view by Connection, the Impulse will be under one of the USB Root Hubs as:

USB Serial Converter / USB Serial Port (COMx).

Note: If Device Manager only lists the USB Serial Converter but not the COM port it could be that the Virtual COM Port driver is not enabled. Open USB Serial Converter Properties and go to Advanced. Check the Load VCP box if it is not already checked and press OK. Then the COM port should show up.

Note: You can changed the COM port number assigned by Windows in Device Manager. Open the Properties for the USB Serial Port (COMx), go to Port Settings and press Advanced. Select the desired COM Port Number from the drop down list box and press OK. To get the device list to show the new COM port number perform a Scan for hardware changes.

Note: If Device Manager says that a COM port number is in use, it may be from another USB device that is no longer being used. You can click through the error message and force it to the number you want.

Note: If you unplug your Impulse, you can still see it in Device Manager by selecting View / Show hidden devices. It will be shown grayed out.

ADVANCED USERS

Advanced users can get more information about the FT232R from the FTDI web site: www.ftdichip.com. You can get new software drivers, application notes, and USB utilities. You can learn how to view your USB connections and load and/or delete all FTDI drivers from your PC. You can get drivers for other operating systems. You can even use their D2XX direct interface API to include in your own custom interface programs if you don't want to use a COM port.

COM PORT SETTINGS

Settings for the COM port should be made by the program that opens and uses the COM port such as:

- Terminal emulation program (HyperTerminal, Tera Term or other)
- Your custom Impulse controller program
- Fluke Ansur Test Automation System program

The settings in Device Manager are usually irrelevant because they are overridden by the controlling program.

The COM port should be set to:

- 115,200 baud
- No parity
- 8 data bits
- 1 stop bit
- Hardware handshaking should be turned on.

COMMAND PROTOCOL

Commands may be sent in upper or lower case.

Commands must be terminated by a Carriage Return (CR) (0x0D) and/or a Line Feed (LF) (0x0A).

Some commands require one or more parameters to be sent with them. Where a command needs parameters, the command is followed by an equal sign and the parameters. Multiple parameters are separated by commas.

Space (SP) (0x20) characters are ignored. The Backspace (BS) (0x08) character removes the previously transmitted character from the command. The Escape (ESC) (0x1B) character erases all previously transmitted characters.

COMMAND RESPONSES

After receiving a command, the Impulse will not store or respond to additional received characters until it has executed the command and responded to it.

The Impulse always responds to a command after it has executed it, by returning a response, terminated by a Carriage Return (0x0D) and a Line Feed (0x0A).

The standard command response is "*", unless other data is to be returned. "*" indicates that the command was understood and executed.

Ver 2.4

A few commands remain active after returning an initial response, as described below.

Incorrect commands return the following error codes:

Code	Description
!	Command empty, no characters
!00	No commands allowed now
!01	Unknown command
!02	Illegal command, not allowed in current mode
!03	Illegal parameter
!04	Receive error
!05	General failure
!06	Option not installed, such as Pacer command sent to Impulse 6000D
!20	Defib data not available
!21	Gas gauge bad read
!24	Data corrupted

LOCAL CONTROL

The Impulse powers up initially under local control by user keys. Then, the only legal command is REMOTE that brings Impulse to remote control.

REMOTE CONTROL

In remote control, Impulse accepts commands and executes them. Some commands set Impulse into special modes. Some commands are only legal in certain modes. The modes are listed in the table:

Mode Mnemonic	Description
MAIN	Main remote mode
DEFIB	Measure defib pulses
PAPULSE	Measure pacer pulse parameters
PASENSE	Test pacer sensitivity
PAREFRACT	Measure pacer refractory periods
ECG	Simulate ECG waves.
ECGPACED	Simulate ECG interactively with pacer.
ECGPERF	Simulate ECG waves for performance testing.
ECGNOISE	Simulate noise on ECG.

The LOCAL command brings Impulse back to local control.

PACER COMMANDS (IMPULSE 7000DP ONLY)

Commands for pacer functions only work with the Impulse 7000DP. The Impulse 6000D will respond to them with the !06 error code.

Ver 2.4

COMMAND SPECIFICATIONS

Note: Unless specified otherwise, commands return *.

GENERAL COMMANDS

REMOTE

Modes: Local control.

Description: Goes to remote control MAIN mode.

LOCAL

Modes: All.

Description: Exits remote control and returns to local control.

IDENT

Modes: All.

Description: Asks for the identification with option and software version number.

Returns: The identification: TBD.

SN

Mode: All.

Description: Retrieves serial number.

Returns: Serial number.

VER

Modes: All.

Description: Asks for the software version number.

Returns: The software version: 3 digits with decimal point, format n.nn.

MODE=MODE

Modes: MAIN.

MODE: The mode to go to: The mnemonic of the mode.

Description: Go to the designated mode.

OMODE

Modes: All.

Description: Queries the mode.

Returns: The current mode mnemonic.

QSET

Modes: All.

Description: Queries the current settings.

Returns: Current settings, separated by commas; settings depend on mode as applicable:

first the mode mnemonic, then settings:

 ${\tt DEFIB, LD=load, WV=wave, CV=wave\ conversion, EA=ECG\ amplitude}$

PAPULSE, BR=brand, IN=input, LD=load

PASENSE, BR=brand, IN=input, LD=load, SH=shape, WD=width, PL=polarity,

EA=ECG amplitude

PAREFRACT, BR=brand, IN=input, LD=load

ECG, WV=wave, EA=ECG amplitude

ECGPACED, WV=wave, TH=threshold, EA=ECG amplitude

ECGPERF, SH=shape, FQ=frequency, EA=ECG amplitude [regular performance wave] ECGPERF, SH=shape, WD=width, RT=rate, EA=ECG amplitude [for R wave pulse]

FLUKE Biomedical Impulse 6000D/7000DP Remote Communications Interface Ver 2.4

ECGNOISE, NF=noise frequency, EW=ECG wave on/off, NA=noise amplitude

EXIT

Modes: All.

Description: Turns off measurement and ECG. Exits the current mode and goes to MAIN

mode.

GLOBAL SETUP COMMANDS

DEFLOAD=LOAD

Model: Impulse 7000DP only

Special: This is a special command only to be used when the Impulse 7000DP is

connected to an Impulse 7010 Defibrillator Selectable Load Accessory. It

tells the 7000DP how to compute the defib pulse parameters.

After disconnecting the 7010, be sure to set the load back to 50 ohms or cycle

power. The 7000DP powers up to 50 ohms.

Mode: All.

LOAD: The defib load in ohms: 3 digits: 025 to 200 by 025.

Description: Sets the load to be used for defib tests when connected to the Impulse 7010

Defibrillator Selectable Load Accessory.

PAINPUT=INPUT

Model: Impulse 7000DP only

Mode: All modes

INPUT: The pacer input: **DEFIB** or **PACER**. **Description:** Sets the input to be used for pacer tests.

PALOAD=LOAD

Model: Impulse 7000DP only

Mode: All.

LOAD: The pacer load in ohms: 4 digits: 0050 to 1500 by 0050.

Description: Sets the load to be used for pacer tests. Connects that load to the pacer jacks.

PABRAND=BRAND

Model: Impulse 7000DP only

Modes: All.

BRAND: The pacer brand: NONE, PHYSIO, PHILIPS, ZOLL, CARDIAC, MRL,

SCHILLER, or **MDE**. (Also **MEDTRONIC** for legacy, same algorithm as **PHYSIO**).

Description: Sets the pacer brand algorithm to be used for pacer tests.

ECGAMPL=AMPLITUDE

Modes: All modes except **PASENSE** and **ECGNOISE**.

AMPLITUDE: The amplitude in mV: 3 digits with decimal point: 0.05 to 5.00.

Description: Sets the ECG wave amplitude for all ECG waves except for Pacer Sensitivity and

ECG Noise waves.

ECGREF=LEAD

Modes: All.

LEAD: The ECG reference lead II or I. **Description:** Sets the ECG reference lead.

DEFIB MODE COMMANDS

DCONVERT=WAVE

Mode: DEFIB.

WAVE: The defib post-shock conversion wave:

CONVERT Convert to normal sinus 60 bpm.

NOCONVERT No change to wave. **ASYSTOLE** Change to asystole wave.

SYNCCONVERT Convert to normal sinus 60 bpm only if sync time within

range of -120 to +380 ms, otherwise change to asystole.

Description: Sets the defib post-shock conversion wave.

DAFIB=GRAN

Mode: DEFIB.

GRAN: The afib granularity: COARSE or FINE.

Description: Runs the afib wave.

DVFIB=GRAN

Mode: DEFIB.

GRAN: The vfib granularity: **COARSE** or **FINE**.

Description: Runs the vfib wave.

DVFIB2=GRAN

Mode: DEFIB.

GRAN: The vfib granularity: **COARSE** or **FINE**.

Description: Runs the vfib2 wave. Vfib2 is a faster version of vfib.

DMONOVTACH=RATE

Mode: DEFIB.

The mono vtach rate in BPM: 3 digits: 120 to 300. Description: Runs the mono vtach wave at the specified rate.

DPOLYVTACH=CODE

Mode: DEFIB.

CODE: The poly vtach code: 1 digit: 1 to 5.

Description: Runs the poly vtach wave of the specified code.

DNSR=RATE

Mode: DEFIB.

RATE: The normal sinus rate in BPM: 3 digits: 150 to 300. **Description:** Runs the normal sinus wave at the specified rate.

DASYSTOLE

Mode: DEFIB.

Description: Runs the asystole wave.

DREADY

Mode: DEFIB.

Description: Ready the measurement system to wait for and measure a defib pulse.

Returns: * (Crlf). Then waits for the defib pulse.

After detecting and processing the defib pulse:

Returns the defib data in numeric fields separated by commas:

```
Type of pulse:
              1 for Monophasic
              2 for Bi-Phasic
              3 for Pulsed Bi-Phasic
For type 1 pulse, the remaining fields are:
       Energy (J):
                             XXX.X
       Peak Voltage (V):
                             XXXX
       Peak Current (A):
                             XXX.X
       Pulse Width 50% (ms): xx.x
       Pulse Width 10% (ms): xx.x
       Sync Time (ms):
                             ±xxx
       ECG Wave now running:
              n for no change.
              c for converted to NSR at 60 bpm.
              A for Asystole.
       Charge Time (s):
                             XXX.X
       Example: 1,123.4,2000,040.2,08.3,12.4,+120,N,012.3
For type 2 pulse, the remaining fields are:
       Energy (J):
                                    XXX.X
       Phase 1 Peak
                       Voltage (V):
                                    XXXX
       Phase 1 Average Voltage (V):
                                    XXXX
       Phase 1 Peak
                       Current (A):
                                    XXX.X
       Phase 1 Average Current (A): xxx.x
       Phase 1 Pulse Width (ms):
                                    XX.X
       Phase 2 Peak
                       Voltage (V):
                                    XXXX
       Phase 2 Average Voltage (V):
                                    XXXX
       Phase 2 Peak
                       Current (A):
                                    XXX.X
       Phase 2 Average Current (A):
                                    XXX.X
       Phase 2 Pulse Width (ms):
                                    XX.X
       Inter-Phase Delay (ms):
                                    XX.X
       Tilt (%):
                                    XX
       Sync Time (ms):
                                    ±xxx
       ECG Wave now running:
              n for no change.
              c for converted to NSR at 60 bpm.
              A for Asystole.
       Charge Time (s):
                                    XXX.X
Example: 2,123.4,2000,1453,040.2,033.1,10.3,
                  1256,0967,032.2,018.1,09.2,
                   02.3,12,+120,N,012.3
For type 3 pulses, the remaining fields are:
       Energy (J):
                                    XXX.X
       Phase 1 Peak
                       Voltage (V):
                                    XXXX
       Phase 1 Average Voltage (V):
                                    XXXX
       Phase 1 Peak
                       Current (A):
                                    XXX.X
       Phase 1 Average Current (A): xxx.x
       Phase 1 Pulse Width (ms):
                                    XX.X
```

Phase 2 Peak Voltage (V): XXXX Phase 2 Average Voltage (V): XXXX Phase 2 Peak Current (A): XXX.X Phase 2 Average Current (A): XXX.X Phase 2 Pulse Width (ms): XX.X Inter-Phase Delay (ms): XX.X Tilt (%): XX Frequency (Hz): XXXX Duty Cycle (%): XX Sync Time (ms): **±**XXX

ECG Wave now running: **n** for no change.

c for converted to NSR at 60 bpm.

A for Asystole.

Charge Time (s): xxx.x

Example: 3,123.4,2000,1453,040.2,033.1,10.3, 1256,0967,032.2,018.1,09.2, 02.3,12,4023,41,+120,N,012.3

Exit: If no defib pulse comes, will guit after receiving ESC character, then sends Crlf

and quits.

DWAVEDATA

Mode: DEFIB.

Description: Requests the defib pulse data from the last measured defib pulse.

Returns: The defib pulse data: 2,500 signed current readings, in Amps, 20 µs apart,

formatted: **±xxx.x**, separated by commas. With a CRLF after every 10 readings.

Example: +001.2,+002.3,-043.2,+100.0, ...

ErrorMessage: Returns !20 if no defib pulse data available.

PACER PULSE COMMANDS

PAREADY

Model: Impulse 7000DP only

Modes: PAPULSE.

Description: Ready the measurement system to wait for and measure pacer pulses

continuously.

Returns: *. Then waits for pacer pulses. After processing each pacer pulse, returns the

pacer data in numeric fields separated by commas:

Rate (PPM): xxx.x
Pulse Width (ms): xxx.xx
Energy (uJ): xxxxxx
Amplitude (mA): +xxx.xx

Example: 120.4,021.63,0146343,+118.62

It takes 2 pulses to calculate a rate. For the 1st pulse, the rate will be returned as

000.0.

Continues sending pacer data for every pacer pulse detected.

Note: For version 2.02 and after, amplitude is always positive. We continue

to use the '+' sign for legacy interface code.

Exit: Will quit after receiving **ESC** character, then sends **Crlf** and quits.

PACER SENSITIVITY COMMANDS

PASRWAVE=WAVE, WIDTH, POLARITY

Model: Impulse 7000DP only

Modes Pasense.

WAVE: The wave shape: 3 characters: **FLT** for flat (off), **SQR** for square, **TRI** for triangle,

SIN for sine.

WIDTH: The width in ms: 3 digits: 001 to 300.

POLARITY: The polarity: 0 for positive, 1 for negative.

Description: Runs the specified pacer sensitivity test pulse.

PASAMPL=AMPLITUDE

Model: Impulse 7000DP only

Mode: PASENSE.

AMPLITUDE: The amplitude in mV: 3 digits with decimal point: 0.05 to 5.00.

Description: Sets the pacer sensitivity wave amplitude.

PASAUTO

Model: Impulse 7000DP only

Mode: PASENSE.

Description: Runs the sensitivity test automatically. This test takes several seconds

depending on the rate. The test interacts with pacer pulses to determine the

sensitivity threshold amplitude for the sensitivity wave that is running.

Returns: * (Crlf). Then waits for pacer pulses and returns intermediate and final test

data for amplitude. Returns data:

R=XXX Pacing rate measured from the initial 3 pulses, only once.

A~X.XX Intermediated amplitude, every pacer pulse.

A=X.XX Final amplitude, only once.

Exit: This test will exit before completion if it receives any alphabetic character. Then it

returns *.

PACER REFRACTORY COMMANDS

PARAUTO

Model: Impulse 7000DP only

Mode: PAREFRACT.

Description: Runs the refractory period test automatically. This test takes several seconds

depending on the pacing rate. The test interacts with pacer pulses to determine the Pulse Refractory Period (PRP), then the Sense Refractory Period (SRP).

Returns: * (Crlf). Then waits for pacer pulses and returns intermediate and final test

data for pacing rate, PRP and SRP. Returns data after every pacer pulse:

R=XXX Pacing rate, only once, after 2 pulses.

P~XXX Intermediate PRP, every pacer pulse after rate.

P=XXX Final PRP, only once.

s~xxx Intermediate SRP, every pacer pulse after final PRP.

s=xxx Final SRP, only once.

Exit: This test will exit before completion if it receives any alphabetic character. Then it

returns *.

ECG NOISE COMMANDS

NOISE=FREQUENCY, ECG

Modes: ECGNOISE.

FREQUENCY: The noise frequency in Hz: 50 or 60.

ECG: Flag to turn ECG NSR60 wave on or off, added to the noise: **T** (on) or **F** (off).

Description: Runs the specified noise wave.

NOISEAMPL=AMPLITUDE

Modes: ECGNOISE.

AMPLITUDE: The noise amplitude in mV: 3 digits with decimal point: 00.0 to 10.0.

Description: Sets the noise amplitude.

ECG COMMANDS

ATRPACE=WIDTH, AMPLITUDE

Mode: ECG.

WIDTH: The width in ms: 2 digits with decimal: 0.1, 0.2, 0.5, 1.0, or 2.0.

AMPLITUDE: The amplitude in mV: 3 digits with sign: ±000 to ±700 (for zero, 0000 also valid).

Description: Sets the width and amplitude for the simulated TV atrial pacer pulse.

VENTPACE=WIDTH, AMPLITUDE

Mode: ECG.

WIDTH: The width in ms: 2 digits with decimal: 0.1, 0.2, 0.5, 1.0, or 2.0.

AMPLITUDE: The amplitude in mV: 3 digits with sign: ±000 to ±700 (for zero, 0000 also valid).

Description: Sets the width and amplitude for the simulated TV ventricular pacer pulse.

NSR=RATE

Mode: ECG.

RATE: The normal sinus rate in BPM: 3 digits: 030 to 360. **Description:** Runs the normal sinus wave at the specified rate.

AFIB=GRAN

Mode: ECG.

GRAN: The afib granularity: COARSE or FINE.

Description: Runs the afib wave.

VFIB=GRAN

Mode: ECG.

GRAN: The vfib granularity: **COARSE** or **FINE**.

Description: Runs the vfib wave.

VFIB2=GRAN

Mode: ECG.

GRAN: The vfib granularity: **COARSE** or **FINE**.

Description: Runs the vfib2 wave. Vfib2 is a faster version of vfib.

MONOVTACH=RATE

Mode: ECG.

RATE: The vtach rate in BPM: 3 digits: 120 to 300. **Description:** Runs the vtach wave at the specified rate.

FLUKE Biomedical Impulse 6000D/7000DP Remote Communications Interface Ver 2.4

POLYVTACH=CODE

Mode: ECG.

CODE: The poly vtach code: 1 digits: 1 to 5.

Description: Runs the poly vtach wave of the specified code.

SPVWAVE=WAVE

Mode: ECG.

WAVE: The supraventricular wave to run:

AFL Atrial Flutter
SNA Sinus Arrhythmia
MBT Missed Beat
ATC ATach

PAT Paroxysmal ATach
NOD Nodal Rhythm
SVT Supra VTach

Description: Runs the supraventricular ECG wave.

PREWAVE=WAVE

Mode: ECG.

WAVE: The premature wave to run:

Atrial PAC PAC Nodal PNC PNC PVC1 **PVC1 Left Vent** PVC1 LV Early PVC1E PVC1 LV R on T PVC1R **PVC2** Right Vent PVC2 PVC2 RV Early PVC2E PVC2R PVC2 RV R on T Multifocal PVCs MF

Description: Runs the premature ECG wave.

VNTWAVE=WAVE

Mode: ECG.

WAVE: The ventricular wave to run:

PVCs 6/min PVC6M PVCs 12/min PVC12M PVCs 24/min PVC24M Freq Multifocal **FMF** Trigeminy TRIG Bigeminy BIG Pair PVCs PAIR Run 5 PVCs RUN5 RUN11 Run 11 PVCs

ASYS Asystole

Description: Runs the ventricular ECG wave.

FLUKE Biomedical Impulse 6000D/7000DP Remote Communications Interface Ver 2.4

CNDWAVE=WAVE

Mode: ECG.

WAVE: The conduction wave to run:

1DB 1º Block 2DB1 2º Block Type I 2DB2 2º Block Type II

3DB 3º Block RBBB RBBB LBBB LBBB

Description: Runs the conduction ECG wave.

TVPWAVE=WAVE

Mode: ECG.

WAVE: The TV paced wave to run:

ATR Atrial 80 BPM
ASY Async 75 BPM
DFS Demand Freq Sinus
DOS Demand Occ Sinus
AVS AV Sequential
NCP Non-Capture
NFN Non-Function

Description: Runs the TV paced ECG wave.

ECG PACED COMMANDS

EPATHRESH=THRESHOLD

Model: Impulse 7000DP only

Mode: ECGPACED.

THRESHOLD: Pacer response threshold mA: 3 digits: 000 to 250 (000 turns off threshold

check and allows all pacer pulses to trigger).

Description: Sets the threshold of pacer amplitude to trigger pacer response wave for pacer

interactive ecg waves.

EPAWAVE=WAVE

Model: Impulse 7000DP only

Mode: ECGPACED.

WAVE: The ECG paced wave to run:

ASY Asynchronous NCP Non-Capture NFN Non-Function

Description: Runs the wave.

EPADEMAND=RATE

Model: Impulse 7000DP only

Mode: ECGPACED.

RATE: The normal sinus rate in BPM for the Demand wave in this mode: 3 digits:

030 to 360.

Description: Runs the Demand wave at the rate.

ECG PERFORMANCE COMMANDS

EPFWAVE=WAVE, FREQUENCY

Modes ECGPERF.

WAVE: The wave shape: 3 characters: FLT for flat (off), SQR for square, TRI for triangle,

SIN for sine.

FREQUENCY: The frequency in Hz: 3 digits with no decimal point 001 to 200; or 4 digits with

decimal point 0.050 to 9.999.

Description: Runs the specified performance wave.

EPFRWAVE=WAVE, WIDTH, RATE

Modes ECGPERF.

WAVE: The wave shape: 3 characters: **FLT** for flat (off), **SQR** for square, **TRI** for triangle,

SIN for sine.

WIDTH: The width in ms: 3 digits: 001 to 300.

RATE: The rate in BPM: 3 digits: 030 to 300.

Pescription: Runs the specified performance pulse.