

05-437

PRIMALERT® 35 Area Radiation Monitor

Technical Data



The PRIMALERT 35 Area Radiation Monitor contains an energy-compensated GM detector and has six range indicators (1, 2, 4, 8, 16, and 32 mR/hr) that can clearly display an increase or decrease in radiation levels. The visible and audible alarms can be set at any of the six levels by a front-panel, screwdriver-adjustable control. When each preset radiation level is exceeded, personnel are alerted by bright flashing red lights (visible over a 180° field) and a loud intermittent audio signal. The alarms stop automatically when the radiation level falls below each of the preset values. This permits instant radiation-level recognition not readily distinguishable on meter-type instruments.

Fail-safe operation is assured by a light which continuously indicates background radiation and provides visual proof that the unit is functioning. The monitor will not jam or show false readings in high radiation fields. A mounting bracket and a 110 V ac adapter/power converter are also included.

Key features

- Provides continuous visual indication of radiation levels and produces audible and/or visual alarms at any of six programmable radiation levels
- Assures reliable, continuous monitoring wherever radioactive materials are present
- Displays the radiation level in bright color-coded lights
- Optional Primalarm Remote Alarm, which functions up to 100 feet from monitor

Specifications

Power requirements	105 V to 125 V/60 Hz/8 W
Dimensions (WxHxT)	9 cm x 15 cm x 4 cm (3.5 in x 6 in x 1.5 in)
Weight	0.9 kg (2 lb)

Ordering information

Model numbers/descriptions

05-437 Primalert 35 Area Radiation Monitor

Standard accessories

AC adapters, specify with order

14-314 110 V ac 12 V dc 500 mA (USA, Japan)

14-400 230 V ac 12 V dc 500 mA (Europe)

14-436 230 V ac 12 V dc 580 mA (Australia)

Optional accessories

62-103 Check Source, ¹³⁷Cs, 10 µCi. Flat disc, 1 in diameter

About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-6 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 certified and our products are:

- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required
- NRC Compliant, where required

Fluke Biomedical.

Better products. More choices. One company.

Fluke Biomedical

6045 Cochran Road
Cleveland, OH 44139-3303 U.S.A.

Fluke Biomedical Europe

Science Park Eindhoven 5110
5692EC Son, The Netherlands

For more information, contact us:

In the U.S.A. (800) 850-4608 or
Fax (440) 349-2307
In Europe/M-East/Africa +31 40 267 5435 or
Fax +31 40 267 5436
From other countries +1 (440) 248-9300 or
Fax +1 (440) 349-2307
Email: sales@flukebiomedical.com
Web access: www.flukebiomedical.com

©2010-2013 Fluke Biomedical. Specifications subject to change without notice. All OEM company trademarks are implied. Printed in U.S.A. 12/2013 3670898C_EN

Modification of this document is not permitted without written permission from Fluke Corporation.