

## Biomedical

# 190M

# Medical ScopeMeter® Portable Oscilloscope

## **Technical Data**



# The 190M: a new generation of medical oscilloscope

The 190M Medical ScopeMeter portable oscilloscope is a high-performance test tool built upon the legacy of Fluke and Fluke Biomedical oscilloscopes in partnership with real customers like you. The 190M is available with choice of two or four channels and offers an unprecedented level of performance, ruggedness, and portability. With the combined power of a high-performance oscilloscope, multimeter and paperless recorder in an easy-to-use test tool, the 190M is the one test tool you can rely on to tackle just about any troubleshooting task in the field.

To minimize downtime and repair costs, you need to get to the root cause of problems as quickly as possible. The 190M offers a number of unique features to help you quickly set up the scope and diagnose difficult problems like intermittent events, signal fluctuations or drift.

Extend your arsenal of troubleshooting capabilities with the new Fluke Biomedical 190M Medical ScopeMeter portable oscilloscope, designed to meet the demands of field service professionals.

## **Key features**

- Two or four electrically-isolated inputs
- Fast sampling rate, up to 2.5 GS/s on two channels simultaneously with up to 400 ps resolution
- Deep memory: 10,000 samples per channel waveform capture so you can zoom in on the details (scope mode)
- Dedicated 5000 count digital multimeter in two-channel model
- Quad meter measurements via scope BNC inputs in four channel model
- Connect-and-View™ triggering for intelligent, automatic triggering on fast, slow and even complex signals
- Frequency spectrum using FFT-analysis
- High-resolution, non-interlaced video
- Smart averaging
- ScopeRecord roll mode gives 30,000 points per input channel and capture waveform sample data for up to 48 hours

- TrendPlot, trend measurement readings for up to 22 days
- Advanced automatic measurements, power (Vpwm, VA, W, PF) and time (mAs, V/s, w/s)
- Two USB ports make it easy to transfer data to a PC and store unlimited waveforms, screen captures and instrument setups on USB memory devices
- New high-perfomance Li-ion battery technology delivers the longest battery life on the market
- Charge spare battery using optional external battery charger
- Easy-access battery door for quick swaps in the field
- Security slot locks down oscilloscope with Kensington lock while unattended
- Environmentally tested to meet IP-51 and withstand 3 g vibration or 30 g shock



## **Technical specifications**

|                               | 190M-2  | 190M-4                                     |
|-------------------------------|---|--|
| Oscilloscope modes            | <u>'</u>  |  |
| Vertical deflection           |   |  |
| Number of channels            | 2   | 4  |
| Bandwidth                     | 200 MHz   | -  |
| Rise time                     | 1.7 ns  |  |
| Number of scope inputs        | 2 input channels plus external trigger  | 4 input channels                           |
| Channel architecture          | All inputs fully insulated from each other  |  |
|                               | activated in any combination  |  |
| Input coupling                | AC or DC, with ground level indicator   |  |
| Input sensitivity             | 2 mV/div to 100 V/div, plus variable atte   | nuation                                    |
| Bandwidth limiter             | User selectable: 20 kHz, 20 MHz or full k   | pandwidth                                  |
| Normal/invert/variable        | On each input channel, switched separa  | tely                                       |
| Extended offset               | Not avalable currently  | -  |
| Input voltage                 | CAT III 1000 V/CAT IV 600 V rated, see g  | eneral specifications for further details  |
| Vertical resolution           | 8 bit   | -  |
| Accuracy                      | ± (2.1 % of reading + 0.04 x range/div)   | @ 5 mV/div to 100 V/div                    |
| Input impedance               | $1 \text{ M}\Omega \pm 1 \text{ %/14 pF} \pm 2 \text{ pF}$  |  |
| Horizontal                    |   |  |
| Maximum real-time sample rate | 2.5 GS/s (2ch)  | 2.5 GS/s (2ch) 1.25 GS/s (4ch)             |
| (sampled simultaneously)      |   |  |
| Record length                 | Up to 10,000 samples per channel  |  |
| Time base range               | 2 ns/div to 4 s/div   |  |
| _                             | Time base in a 1-2-4-sequence Slower time/division settings using   |  |
|                               | ScopeRecord™ roll mode (see recorder mode)  |  |
| Maximum record length         | 10,000 samples per channel in scope mode  |  |
|                               | 30,000 points per channel in ScopeRecord™ roll mode (see recorder mode)   |  |
| Timing accuracy               | ± (0.01 % of reading + 1 pixel)   |  |
| Glitch capture                | 8 ns peak detect on each channel (using real time sampling and data   |  |
| _                             | compression, at any timebase setting)   |  |
| Display and acquisition       |   |  |
| Display                       | 153 mm (6 in) full-color LCD with LED backlight   |  |
| Display modes                 | Any combination of channels; average on/off; replay   |  |
| Visible screen width          | 12 divisions horizontally in scope mode   |  |
| Digital persistence modes     | Off/short/medium/long/infinite and envelope mode  |  |
| Waveform mathematics          | A + B, A - B, A x B, all with user-selectable scaling of resultant;   |  |
|                               | A versus B (X-Y- mode); frequency spectrum using FFT analysis   |  |
| Acquisition modes             | Normal, averaged, auto, single shot, ScopeRecord™ roll, glitch capture,   |  |
| Trigger and dolor-            | waveform compare with automatic pass/fail testing; replay   |  |
| Trigger and delay Source      | Input A, B or external (via meter input)  | Input A,B,C or D                           |
| Modes                         |   | -  |
| INIOGES                       | Automatic Connect-and-View <sup>TM</sup> , free run, single shot, edge, delay, dual slope, video, video line, selectable pulsewidth (channel A only), N-cycle |  |
| Connect-and-View™             |   | gnizes signal patterns, automatically sets |
|                               | up and continuously adjusts triggering, time base and amplitude   |  |
|                               | Automatically displays stable waveforms   |  |
|                               | motor drive and control signals can be s  | witched off if preferred                   |



## ■ Biomedical

|  | 190M-2 190M-4   |  |
|--|---|--|
| Video triggering (on ch. A)              | NTSC, PAL, PAL+, SECAM; includes field 1, field 2 and line select   |  |
| High-res, non-interlaced video           | Non-interlaced video with line-select, for line frequencies in the range 14 kHz up to 65 kHz  |  |
| Pulse width triggering (on channel A)    | Pulse width qualified by time allows for triggering $<$ t, $>$ t, $=$ t, $\neq$ t, where t is selectable in minimum steps of 0.01 div or 50 ns  |  |
| Time delay                               | 1 Full screen of pre-trigger view or up to 100 screens (= 1,200 divisions) of post-trigger delay  |  |
| Dual slope triggering                    | Triggers on both rising and falling edges alike   |  |
| N-cycle triggering                       | Triggers on Nth occurrence of a trigger event; N to be set in the range 2 to 99   |  |
| <b>Automatic capture of 100 screens</b>  |   |  |
| When an anomaly is seen, the replay bu   | ent always memorizes the last 100 screens—no specific user setup required.<br>tton can be pressed to review the full sequence of screen events over and over.<br>es or intermittent anomalies and will operate in baby-sit mode capturing 100 |  |
| Replay                                   | Manual or continuous replay. Displays the captured 100 screens as a live animation or under manual control. Each screen has date and time-stamp   |  |
|  | Two sets of 100 screens each can be saved internally for later recall and analysis Direct storage of additional sets on external flash memory drive through USB host port   |  |
| Fast Fourier Transform (FFT) frequenc    |   |  |
| Shows frequency content of oscilloscope  |   |  |
| Window                                   | Automatic, hamming, hanning or none   |  |
| Automatic window                         | Digitally re-samples acquired waveform to obtain optimum frequency resolution in FFT resultant  |  |
| Vertical scale                           | Linear/logarithmic (in volts or amps)   |  |
| Frequency axis                           | Logarithmic frequency range automatically set as a function of timebase range of oscilloscope   |  |
| Waveform compare and pass/fail testi     | ng  |  |
| Waveform Compare                         | Provides storage and display of a reference waveform for visual comparison with newly acquired waveforms. Reference is derived from an acquired waveform and can be modified in the oscilloscope or externally using FlukeView Software.      |  |
| Pass/Fail Testing                        | In waveform compare mode, the oscilloscope can be set to store only matching (pass) or only non-matching (fail) acquired waveforms in the replay memory bank for further analysis   |  |
| <b>Automatic scope measurements</b>      |   |  |
|  | eak min, Vpeak to peak, A ac, A dc, A ac + dc, frequency (in Hz), rise time phase (between any 2 inputs), pulse width (pos./neg.), duty cycle (pos./neg.), pan), dBV, dBm into 50 $\Omega$ and 600 $\Omega$                                   |  |
| Advanced power and motor drive functions | V/Hz Ratio (190M-2 only), Power Factor (PF), watts, VA, VA reactive, VPWMac and VPWM (ac + dc) for measurement on pulse width modulated motordrives and frequency inverters   |  |
| Advanced functions                       | mA×s (Current-over-time, between cursors); V×s (voltage over time, between cursors); W×s (energy, between cursors)  |  |
| Cursor measurements                      |   |  |
| Source                                   | On any input waveform or on mathematical resultant waveform (Excluding X-Y-mode)  |  |
| Dual horizontal lines                    | Voltage at cursor 1 and at cursor 2, voltage between cursors  |  |
| Dual vertical lines                      | Time between cursors, 1/T between cursors (in Hz), voltage between markers, rise time with markers, fall time with markers; Vrms between cursors, watts between cursors   |  |



## Biomedical

|   | 190M-2  | 190M-4   |  |  |
|---|---|--|--|--|
| Single vertical line  | Min/max and average voltage at cursor position; frequency and rms-value of individual frequency component in the FFT resultant                        |  |  |  |
| ZOOM  | Ranges from full record overview to zoom-in up to sample level at any record length   |  |  |  |
| <b>Meter Modes</b>  |   |  |  |  |
| Meter inputs  | Via 4 mm banana inputs, fully isolated from scope inputs and scope ground   | Via BNC scope inputs   |  |  |
| Number of readings  | One at a time   | Up to 4 simultaneously   |  |  |
| Maximum resolution  | 5,000 counts  | 99 counts  |  |  |
| Input impedance   | $1 \text{ M}\Omega \pm 1 \%/14 \text{ pF} \pm 2 \text{ pF}$   |  |  |  |
| Advanced meter functions  | Auto/manual ranging, relative measurements (Zero reference),<br>TrendPlot™ recording  |  |  |  |
|   | The specified accuracy is valid over the temperature range 18 °C to 28 °C Add 10 % of specified accuracy for each degree C below 18 °C or above 28 °C |  |  |  |
| Voltage   |   |  |  |  |
| Vdc accuracy  | ± (0.5 % + 5 counts)  | ± (0.5 % + 5 counts)   |  |  |
| Vac true rms accuracy 15 Hz to 60 Hz: 60 Hz to 1 kHz: 60 Hz to 20 kHz:    | ± (1 % + 10 counts)<br>± (2.5 % + 15 counts)  | ± (1.5 % + 10 counts)<br>± (2.5 % + 15 counts)                                     |  |  |
| Vac+dc true rms accuracy 15 Hz to 60 Hz: 60 Hz to 1 kHz: 60 Hz to 20 kHz: | ± (1 % + 10 counts)<br>± (2.5 % + 15 counts)  | ± (1.5 % + 10 counts)<br>± (2.5 % + 15 counts)                                     |  |  |
| Voltmeter ranges  | 500 mV, 5 V, 50 V, 500 V, 1,000 V   |  |  |  |
| Resistance  |   |  |  |  |
| Ranges  | 500 Ω, 5 kΩ, 50 kΩ,<br>500 kΩ, 5 MΩ, 30 MΩ  | Feature/function not available for this model                                      |  |  |
| Accuracy  | ± (0.6 % + 5 counts)  |  |  |  |
| Other meter functions   |   |  |  |  |
| Continuity  | Beeper on $<$ 50 $\Omega$ ( $\pm$ 30 $\Omega$ )   | Feature/function not available   |  |  |
| Diode test  | Up to 2.8 V   | for this model   |  |  |
| Current (A)   | A dc, A ac, A ac + dc using an optional current clamp or shunt<br>Scaling factors: 0.1 mV/A, 1 mV/A to 100 V/A and 400 mV/A                           |  |  |  |
| Temperature   | With optional accessories. Scale factors 1 °C/mV or 1 °F/mV   |  |  |  |
| Recorder Modes  |   |  |  |  |
| ScopeRecord™ Roll Mode  | ScopeRecord™ Roll Mode  |  |  |  |
| Dual or multiple input waveform storage                                   | mode, using deep memory   |  |  |  |
| Source and display  | Input A, Input B, Dual All channels sampled simultaneously  | Any combination of inputs, up to four channels All channels sampled simultaneously |  |  |
| Bandwidth   | 20 MHz or 20 kHz, user selectable   |  |  |  |
|   | 20 Mil 01 20 Mil, abel beleetable   |  |  |  |



| Biomedical |  | Bi | om | edi | ica | ı |
|------------|--|----|----|-----|-----|---|
|------------|--|----|----|-----|-----|---|

|  | 190M-2  | 190M-4                                   |  |
|--|---|--|--|
| Memory depth                               | 30,000 data points, each holding min/m  | nax pair of information                  |  |
| Min/max values                             | Min/max values are created at samples that are measured at high sample rate,  |  |  |
|  | ensuring capture and display of glitches  |  |  |
| Recording modes                            | Single sweep, continuous roll, Single sweep, continuous roll  |  |  |
|  | Start-on-trigger (through external),  | Start-on-trigger (through any channel),  |  |
|  | Stop-on-trigger (through external)  Stop-on-trigger (through any channel  |  |  |
| Stop-on-trigger                            | ScopeRecord mode can be stopped by an individual trigger event or by an   |  |  |
|  | interruption of a repetitive trigger signal through any input channel (through  |  |  |
| Transportationals                          | external on 190M-2 model)   |  |  |
| Horizontal scale                           | Time from start, time of day  |  |  |
| Zoom                                       | Ranges from full record overview to zoom in up to sample level, at any record length  |  |  |
| Memory                                     | Two multiple input ScopeRecord wavefo   | rms can be saved internally for later    |  |
|  | recall and analysis. Direct storage on ex-  | ternal flash memory drive through USB    |  |
|  | host port   |  |  |
| ScopeRecord™ Roll mode sample rate         | <u> </u>  |  |  |
| Time base range                            | 5 ms/div to 2 min/div   |  |  |
| Recorded timespan                          | 6 sec to 48 hr  |  |  |
| Time/division in 'view all' mode           | 0.5 s/div to 4 h/div  |  |  |
| Glitch capture                             | 8 ns  |  |  |
| Sample rate                                | 125 MS/s  |  |  |
| Resolution                                 | 200 μsec to 4.8 sec   | 200 µsec to 4.8 sec                      |  |
| Trendplot™ Recording                       |   |  |  |
|  | corder graphically plots, displays and stor   | es results of up to four automatic scope |  |
| neasurements or a DMM-reading over time    |   |  |  |
| Source and display                         | Any combination of scope measurements, made on any of the input channels, or DMM reading (two-channel instruments)  |  |  |
| Memory depth                               | 18,000 Points (sets) per measurement; each recorded sample point contains a minimum, a maximum and an average value, plus a date and time stamp                 |  |  |
| Ranges                                     | Normal view: 5 s/div to 30 min/div  |  |  |
|  | In view-all mode: 5 min/div to 48 hr/div (overview of total record)   |  |  |
| Recorded time span                         | Up to 22 days, with a resolution of 102 seconds   |  |  |
| Recording mode                             | Continuous recording, starting at 5 s/div with automatic record compression   |  |  |
| Measurement speed                          | 3 Automatic measurements per second or more   |  |  |
| Horizontal scale                           | Time from start, time of day  |  |  |
| Zoom                                       | Up to 64x zoom-out for full record overview, up to 10x zoom-in for maximum detail   |  |  |
| Memory                                     | Two multiple input TrendPlot records can be saved internally for later recall and analysis. Direct storage on external flash memory drive through USB host port |  |  |
| <b>Cursor measurements: all recorder m</b> |   |  |  |
| Source                                     | Any waveform trace in any waveform display mode   |  |  |
|  | (Scope, ScopeRecord or TrendPlot)   |  |  |
| Dual vertical lines                        | Cursors may be used to identify min, max or average value of any datapoint in a   |  |  |
|  | record, with time between cursors, time from start or absolute time   |  |  |
|  |   |  |  |



■ Biomedical

|   | 190M-2 190M-4  |  |
|---|--|--|
| Conorol emogifications  | 190M-4   |  |
| General specifications  |  |  |
| Input voltage range   | GAM III 1000 II/GAM III 000 II   |  |
| Rated maximum floating voltage  | CAT III 1000 V/CAT IV 600 V  |  |
| N   | (Maximum voltage between any contact and earth-ground voltage level)   |  |
| Maximum probe voltage   | CAT III 1000 V/CAT IV 600 V  |  |
| Marrian DNC input roltogo   | (Maximum voltage between any contact and earth-ground voltage level)   |  |
| Maximum BNC input voltage   | CAT IV 300 V (Maximum voltage on BNC input directly)   |  |
| Maximum voltage on meter input  | designed banana input connectors)  |  |
|   | designed banana input connectors)  |  |
| Memory save and recall  | 10777  |  |
| Memory locations (internal)   | 15 Waveform memories plus 2 recording memories   |  |
| 15 waveform memory locations  | Stores ScopeTrace waveform data (2 traces each) plus screen-copy plus  |  |
| Maria de la Caración | corresponding setup  |  |
| Two recording memories  | Each may contain:  |  |
|   | • a 100-screen replay sequence, or   |  |
|   | <ul> <li>a ScopeRecord roll-mode recording (two traces), or</li> <li>a TrendPlot recording of up to four measurements</li> </ul> |  |
| External data storage   | On PC, using FlukeView™ Software, or   |  |
| Likernar data storage   | Direct storage on external flash memory drive (maximum 2 GB) through USB   |  |
|   | host port  |  |
| Screencopies  | On PC, using FlukeView™ Software, or   |  |
| boroomoopios  | • Internally (in instrument), which can be copied on to external flash memory  |  |
|   | drive as .BMP-file through USB host port   |  |
| Volatility  | Measurement data is initially stored in RAM, which is maintained by the main   |  |
|   | battery with a 30-seconds back-up when battery is exchanged  |  |
|   | When storing data, this is written in non-volatile flash-ROM   |  |
| Real-time clock   | Provides date and time stamp information for ScopeRecord,  |  |
|   | for 100-screen replay sequences and for TrendPlot recordings   |  |
| Case  |  |  |
| Design  | Rugged, shock-proof with integrated protective holster.  |  |
|   | Handstrap and hangstrap included as standard Kensington lock supported to lock   |  |
|   | down instrument when left unattended   |  |
| Drip and dust proof   | IP 51 according to IEC 529   |  |
| Shock and vibration   | Shock 30 g, vibration (sinusoidal) 3 g according to MIL-PRF-28800F Class 2   |  |
| Display size  | 127 mm x 88 mm (153 mm/6.0 in diagonal) LCD  |  |
| Resolution  | 320 x 240 pixels   |  |
| Contrast and brightness   | User adjustable, temperature compensated   |  |
| Brightness  | 200 cd/m² typical using power adapter, 90 cd/m² typical using battery power  |  |
| Mechanical data   |  |  |
| Size (HxWxD)  | 265 mm x 190 mm x 70 mm (10.4 in x 7.5 in x 2.8 in)  |  |
| Weight (including battery)  | 2.1 kg (4.6 lb) 2.2 kg (4.8 lb)  |  |
| Power   |  |  |
| Line power  | Mains adapter/battery charger BC190 included, version depending on country   |  |
| Battery power   | Rechargeable double capacity Li-Ion battery (included). Battery swappable  |  |
|   | through easily-accessible battery door at the rear of the instrument   |  |
|   |  |  |



|                                       | 190M-2   | 190M-4                                      |  |
|---------------------------------------|--|---|--|
| Battery type (included) and capacity  | BP290; 2400 mAh [BP291 (4800 mAh)  | BP291; 4800 mAh                             |  |
| [+opt. battery]                       | optional]  | Di 201, 4000 limi                           |  |
| Battery charge indicator              | Battery has built-in status indicator for use with external charger, next to battery                                   |   |  |
| Success on angle mandator             | status indicator on instrument screen  | 50 Willi Oxformar offargor, from to battory |  |
| Battery operating time (with          | Up to four hours using BP290   | Up to seven hours using BP291               |  |
| backlight low)                        | (included); up to eight hours using  | (included)                                  |  |
| ,                                     | BP291 (optional)   | ,   |  |
| Battery charging time                 | 2.5 hours using BP290; 5 hours using   | 5 hours BP291                               |  |
|                                       | BP291  |   |  |
| Battery power saving functions        | Auto power-down with adjustable power-down time; auto display off with   |   |  |
|                                       | adjustable power-down time; on-screen  | battery power indicator                     |  |
| Safety                                |  |   |  |
| Compliance                            | EN 61010-1:2001, Pollution Degree 2; CAN/CSA C22.2, No. 61010-1-04, with   |   |  |
|                                       | approval; UL61010B; ANSI/ISA-82.02.01  |   |  |
| Environmental                         |  |   |  |
| Operating temperature                 | 0 °C to +40 °C; +40 °C to +50 °C Exclud  | ling battery                                |  |
| Storage temperature                   | -20 °C to +60 °C   |   |  |
| Humidity                              | 10 °C to +30 °C: 95 % RH Non-condensing  |   |  |
|                                       | 30 °C to +40 °C: 75 % RH Non-condensing  |   |  |
|                                       | 40 °C to +50 °C: 45 % RH Non-condensing  |   |  |
| Maximum operating altitude            | Up to 2,000 m (6666 ft) for CAT IV 600 V, CAT III 1000 V; up to 3,000 m (10,000  |   |  |
|                                       | ft) for CAT III 600 V, CAT II 1000 V   |   |  |
| Maximum storage altitude              | 12 km (40,000 ft)  |   |  |
| Electro-magnetic-compatibility (EMC)  | EN 61326 (2005-12) For emission and immunity   |   |  |
| Interfaces                            | Two USB ports provided. Ports are fully insulated from instrument's floating   |   |  |
|                                       | measurement circuitry.   |   |  |
|                                       | USB-host port directly connects to external flash memory drive (up to 2 GB for   |   |  |
|                                       | storage of waveform data, complete datasets in which data and setup information  |   |  |
|                                       | is included, instrument settings and screen copies.  |   |  |
|                                       | A mini-USB-B is provided which allows for interconnection to PC for remote control and data transfer under PC-control. |   |  |
| Probe calibration output              |  |   |  |
| Tiobe campiation output               | Dedicated probe-cal output with reference contact provided, fully insulated from any measurement input channel         |   |  |
| Warranty                              | Three-years (parts and labor) on main instrument; one-year on accessories  |   |  |
| Included accessories                  | Timee years (parts and taber) on main instrument, one-year on accessories  |   |  |
| Batterey charger/mains adapter        | BC190  |   |  |
| Li-Ion battery pack                   | BP290 (2400 mAh)   | BP291 (4800 mAh)                            |  |
| Voltage probe sets. Each set includes | VPS410 (One red, one blue)   | VPS410 (One red, one grey, one blue,        |  |
| ground lead, hook clip, ground        | vi bito (one rea, one blue)  | one green)                                  |  |
| spring and probe tip insulation       |  | 910011                                      |  |
| sleeve                                |  |   |  |
| Test leads                            | TL175 (One red, one black) with test   | N/A   |  |
|                                       | pins   |   |  |
| Other                                 | Handstrap affixed to instrument; hangstrap (user-selectable for left- or   |   |  |
|                                       | right-hand use); multi-language users manuals on CD-ROM; FlukeView® demo   |   |  |
|                                       | package (with restricted functionality); USB interface cable for PC connectivity                                       |   |  |
|                                       | F 90 ( 10001000 Tallottoliality), 0  |   |  |



#### Biomedical

## **Ordering information**

## Item numbers/descriptions

190M-2 Medical ScopeMeter Portable Oscilloscope

#### **Included accessories:**

VPS410-R Voltage probe set, 10:1, 300 MHz, one set red

VPS410-B Voltage probe set, 10:1, 300 MHz, one

**TL175** TwistGuard<sup>™</sup> safety-designed test leads set (1 red, 1 black)

EBC290 External battery charger for BP290 and BP291

**SW90W** FlukeView Software for Windows (full version)

C290 Hard shell protective carrying case for 190 Series II

BP290 Li-Ion battery pack, 2400 mAh MA190 Medical Accessory Kit (includes 50 ohm BNC feed-through, 50 ohm 10:1 attenuator feed through, 1 ohm current shunt, 50 ohm current shunt, 50 ohm coax cable, female BNC to 4 mm banana adapter, two female to female 4 mm banana plug adapters)

190M-4 Medical ScopeMeter Portable Oscilloscope

#### **Included accessories:**

VPS410-R Voltage probe set, 10:1, 300 MHz, one set red

VPS410-G Voltage probe set, 10:1, 300 MHz, one set grey

VPS410-B Voltage probe set, 10:1, 300 MHz, one

VPS410-V Voltage probe set, 10:1, 300 MHz, one set green

**EBC290** External battery charger for BP290 and

**SW90W** FlukeView Software for Windows (full

**C290** Hard shell protective carrying case for 190 Series II

BP291 Li-Ion battery pack, 4800 mAh MA190 Medical Accessory Kit (includes 50 ohm BNC feed-through, 50 ohm 10:1 attenuator feed through, 1 ohm current shunt, 50 ohm current shunt, 50 ohm coax cable, female BNC to 4 mm banana adapter, two female to female 4 mm banana plug adapters)

### **About Fluke Biomedical**

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance.

Highly credentialed and equipped with a NVLAP Lab Code 200566-6 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service

for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

## Fluke Biomedical Regulatory Commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 and ISO 13485 medical device certified and our products are:

- CE Certified, where required
  NIST Traceable and Calibrated
  UL, CSA, ETL Certified, where required

## Fluke Biomedical.

Better products. More choices. One company.

#### Fluke Biomedical

6045 Cochran Road Cleveland, OH 44139-3303 U.S.A.

## Fluke Biomedical Europe

Science Park Eindhoven 5110 5692EC Son. The Netherlands

### For more information, contact us:

In the U.S.A. (800) 850-4608 or Fax (440) 349-2307 In Europe/M-East/Africa +31 40 267 5435 or Fax +31 40 267 5436 From other countries +1 (440) 248-9300 or Fax +1 (440) 349-2307 Email: sales@flukebiomedical.com Web access: www.flukebiomedical.com

©2012 Fluke Biomedical. Specifications subject to change without notice. All OEM trademarks are implied. Printed in U.S.A. 08/2012 4122175B D-EN-N

Modification of this document is not permitted without written permission from Fluke Corporation.