

Biomedical

PS410 Patient Simulator

Technical Data



The PS410 is a compact high-performance simulator for patient monitor testing.

This handheld device simulates a full range of cardiac rhythms and a wide variety of ECG conditions. It includes pacemaker simulation, 35 arrhythmia selections, and adult and pediatric normal-sinus rhythms.

Small enough to fit in a pocket, the handy PS410 weighs less than a pound and is easy to operate. Technicians simply connect the simulator to the device under test and use the PS410 keypad to enter the code presets. The simulator then transmits the selected preset simulations to the device.

Key features

- Handheld
- 12-lead ECG simulation
- 12 arrhythmia selections
- Universal ECG jacks
- Auto sequencing of performance waveforms
- Battery operated
- PS410/DPM1B Bundle kit with custom carrying case for quick ECG/NIBP patient monitor testing



Specifications

Normal sinus rhythm (12-lead w	vith independent outputs referenced to R	L)	
Normal rate	80 BPM		
Selectable rates	30, 40, 60, 80, 100, 120, 140, 160, 300 BPM	180, 200, 220, 240, 260, 280, and	
Accuracy	± 1 %		
Output impedance	940 Ω between leads		
ECG amplitudes	0.5, 1, and 2 mV		
Amplitude accuracy	± 2 % Lead II		
High-level output	1000x Lead II		
Adult or pediatric ECG waveform			
ECG performance (Lead II)			
Square wave	2, and 0.125 Hz		
Pulse	30, 60, and 120 BPM; 60 ms pulse width		
Sine waves	0.5, 5, 10, 40, 50, and 60 Hz (1 mV amplitude only, lead II)		
Triangle wave	2 Hz		
ST segment analysis			
Elevated or depressed	-0.6 mV to 6 mV in, 0.2 mV steps		
Pacemaker			
Pacemaker rhythm	Demand pacer with	Demand pacer with	
Pacer non-capture* Pacer non-function	occasional sinus	frequent sinus	
		A-V sequential	
Artifact selections			
50 Hz artifact	Muscle artifact	Respiration artifact	
60 Hz artifact	Baseline artifact		
Arrhythmia selections			
PVC1 left ventricular focus*	Trigeminy	Nodal rhythm	
PVC1 early, LV focus*	PVCs at a rate of 6 per minute	Irregular rhythm	
PVC1 R on T, LV focus*	PVCs at a rate of 12 per minute	Atrial flutter	
Pair of PVCs*	PVCs at a rate of 24 per minute	Atrial fibrillation	
Run of 5 PVCs*	Ventricular tachycardia	Atrial fibrillation 1/2 size	
Run of 11 PVCs*	Ventricular fibrillation	Atrial tachycardia	
Multi-focal PVCs*	Ventricular fibrillation 1/2 size	First-degree block	
Frequent multi-focal PVCs*	Supraventricular tachycardia	Second-degree block	
PVC2 right ventricular focus*	Premature atrial contraction*	Third-degree block	
PVC2 early, RV focus*	Premature nodal contraction*	Right-bundle-branch block	
PVC2 R on T, RV focus*	Asystole	Left-bundle-branch block	
Bigeminy	Missed beat*		
*Event occurs once. To repeat, enter the	selection again.		



General specifications			
Operation	Outputs are easily selected or adjusted using top-panel controls and the large two-digit numeric display		
Controls	Display	2-line LCD	
	Control	6 control keys; ON/OFF power switch	
Power	9 V battery/battery eliminator		
Housing	ABS plastic case		
Dimensions	11.3 cm L x 7 cm W x 3.4 cm H (6.1 in L x 3.7 in W x 1.3 in H)		
Weight	0.4 kg (0.9 lb)		
Temperature	Operating	15 °C to 35 °C (59 °F to 95 °F)	
	Storage	0 °C to 50 °C (32 °F to 122 °F)	

Ordering information

Item numbers/descriptions

PS410 PS410 Patient Simulator

PS410/DPM1B Bundle PS410/DPM1B Bundle Kit (includes PS410, DPM1B, all accessories, and a custom carrying case)

Standard accessories

MANUAL PS410 Users Manual (printed)

CD-ROM PS410 Users Manual (CD)

BE-UNVSL-IEC320C14 Battery Eliminator 100 V ac to 240 V ac

9 V Battery



Optional accessories

17024 Universal Banana Adapter

17191 Carrying Case, single pocket

17192 Carrying Case, double pocket

About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-6 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001

- CE Certified, where required
 NIST Traceable and Calibrated
 UL, CSA, ETL Certified, where required

Fluke Biomedical.

Better products. More choices. One company.

Fluke Biomedical

6045 Cochran Road Cleveland, OH 44139-3303 U.S.A.

Fluke Biomedical Europe

Science Park Eindhoven 5110 5692EC Son, The Netherlands

For more information, contact us:

In the U.S.A. (800) 850-4608 or Fax (440) 349-2307 In Europe/M-East/Africa +31 40 267 5435 or Fax +31 40 267 5436 From other countries +1 (440) 248-9300 or Fax +1 (440) 349-2307 Email: sales@flukebiomedical.com Web access: www.flukebiomedical.com

©2006-2010 Fluke Biomedical. Specifications subject to change without notice. Printed in U.S.A. 9/2010 2818121B D-EN-N

Modification of this document is not permitted without written permission from Fluke Corporation.