

**Biomedical** 

# Technical data Impulse 6000D/7000DP Defibrillator/External Pacer Analyzer

The Impulse 6000D Defibrillator Analyzer and Impulse 7000DP Defibrillator/Transcutaneous Pacer Analyzer Test Systems are rugged, portable precision test instruments that ensure proper operation and ultimate performance of critical lifesupport cardiac-resuscitation equipment.

The Impulse 6000D and Impulse 7000DP testcapabilities encompass the spectrum of worldwideestablished pulse shapes, showcase breakthrough AED technology compatibility, and outperform in accuracy and standards. Additionally, the Impulse 7000DP incorporates the tests and the extensive range of test loads and measurement algorithms needed to test external transcutaneous pacemakers.

In conjunction with an Impulse 7000DP, the Impulse 7010 Defibrillator Selectable Load Accessory provides multiple loads of 25  $\Omega$ , 50  $\Omega$ , 75  $\Omega$ , 100  $\Omega$ , 125  $\Omega$ , 150  $\Omega$ , 175  $\Omega$ , and 200  $\Omega$  for defibrillator performance testing. A standard USB interface enables computer control and data transfer, and optional Ansur PC-based automation software increases productivity by outfitting users with an easy-to-use method to standardize testing procedures and capture, print and document data.



### **Key features**

- Impulse 7010 Defibrillator Selectable Load Accessory provides multiple loads of 25  $\Omega$ , 50  $\Omega$ , 75  $\Omega$ , 100  $\Omega$ , 125  $\Omega$ , 150  $\Omega$ , 175  $\Omega$ , and 200  $\Omega$  to comply with IEC 60601-2-4 standard (optional)
- Lown, Edmark, trapezoidal, biphasic and pulsed-biphasic defibrillation technology compatibility
- AED technology compatibility
- First-class measurement accuracy  $\pm 1\%$  of reading + 0.1 J
- Intuitive user interface and backlight, easy-to-read display
- Portable, rugged, easy to carry
- Long-lasting, rechargeable battery
- Pacer brand selections
- Pacer input protected against defibrillator output (7000DP only)
- 10 independent ECG outputs that provide 12 lead combinations for standardized clinical signals
- Flexible heart-rate settings (1 BPM step) facilitate rate meter accuracy and alarm testing
- DSP-based measurements enable future firmware and waveforms upgrade
- Unique integrated posts for secure connections
- Two-year extended warranty (no-cost extended warranty available after first-year calibration at any Fluke Biomedical authorized service center)
- Optional Ansur test automation software to standardize testing procedures, capture waveforms and test results, and print and document test results
- Designed, tested, and built to incomparable Fluke quality standards

FLUKE 🛛

LANDAUER® 📈 💛 RaySafe®



## Specifications

General specifications	
Operating temperature	10 °C to 40 °C (50 °F to 104 °F)
Storage temperature	-20 °C to 60 °C (-4 °F to 140 °F)
Humidity	10 % to 90 % non-condensing
Display	LCD display
Communications	USB device port for computer control
Modes of operation	Manual and remote
Power	Internal rechargeable NiMH battery pack for nine hours (typical) operation after full charge or the battery charger can operate the analyzer and charge the battery simultaneously
Battery charger	100 V to 240 V input, 15 V/1.5 A output. For best performance, the battery charger should be connected to a properly grounded ac receptacle
Enclosure	ABS plastic housing
Dimensions (WxDxH)	32 cm x 24 cm x 13 cm (13 in x 9.5 in x 5 in)
Weight	3.02 kg (6.6 lb, 0.1 oz)
Safety standards	CE: IEC/EN61010-1 2nd Edition; Pollution degree 2 CAN/CSA-C22.2 No 61010-1; UL61010-1 C-Tick: Australian EMC
Electromagnetic compatibility standards (EMC)	European EMC: EN61326-1





## Defibrillator analyzer technical specifications

Energy output measurement		
Compatible defibrillator wave- shapes	Lown, Edmark, trapezoidal, dc biphasic, and ac pulsed biphasic Note: AC pulsed biphasic waveform has not been approved in the United States.	
Autoranged measurement	0.1 J to 600 J	
Accuracy	0.1 J to 360 J: $\pm$ (1 % of reading + 0.1 J) 360 J to 600 J: $\pm$ (1 % of reading + 0.1 J), typical Note: For pulsed biphasic defibrillator, specified accuracy is $\pm$ (1.5 % of reading + 0.3 J) on both ranges.	
Load resistance		
Resistance	50 Ω	
Accuracy	$\pm$ 1 %, non-inductive (< 2 $\mu\text{H})$	
Pulse trigger level	20 V	
Pulse width		
Range	1 ms to 50 ms	
Accuracy	± 0.1 ms	
Voltage		
Range	20 V to 5000 V	
Accuracy	$\pm$ (1 % of reading + 2 V)	
Current		
Range	0.4 A to 100 A	
Accuracy	± (1 % of reading + 0.1 A)	
Tilt (biphasic and pulsed biphasic)		
Range	1% to 99%	
Accuracy	±1 digit	
Interphase delay (biphasic and pulsed biphasic)		
Range	0.1 ms to 9.9 ms	
Accuracy	± 0.1 ms	
Frequency (pulsed biphasic only)		
Range	2000 Hz to 8000 Hz	
Accuracy	±1% of reading	
Duty cycle (pulsed biphasic only)		
Range	1% to 99%	
Accuracy	± 1 digit	
Sample rate	250 kHz (4 μs sample)	
Maximum average power	12 W, equivalent to 10 defib pulses of 360 J every 5 minutes	
Scope output	Autorange: 2000:1, 400:1, and 80:1 depending on range	
Waveform playback	Output: BNC	
	Output impedance: 50 Ω	
	Amplitude accuracy: ± 5 %	



Charge time measurement		
Range	0.1 s to 100 s	
Accuracy	± 0.05 s, typical	
Synchronization test (elective cardioversion)		
Delay time measurement	<ul> <li>Timing window: ECG R-wave peak to the defib pulse peak</li> <li>Range: -120 ms to 380 ms; measures timing from 120 ms prior to the R-wave peak to up to 380 ms following the R-wave peak</li> <li>Resolution: 1 ms</li> <li>Accuracy: ± 1 ms</li> </ul>	
ECG waves	<ul> <li>Normal sinus rhythm (NSR): 10 BPM to 180 BPM in 1 BPM steps</li> <li>Atrial fibrillation: Coarse and fine</li> <li>Monomorphic ventricular tachycardia: 120 BPM to 240 BPM in 5 BPM steps</li> <li>Asystole: Flat line</li> </ul>	
Automated defibrillator test ECG waves		
Normal sinus	10 BPM to 300 BPM in 1 BPM steps	
Ventricular fibrillation	Coarse and fine	
Monomorphic ventricular tachycardia	120 BPM to 300 BPM in 5 BPM steps	
Polymorphic ventricular tachycardia	5 types	
Asystole	Flat line	

## ECG waves

ECG general		
Lead configuration	12-lead simulation; RA, LL, LA, RL, V1-6 with independent outputs Lead to lead impedance: 1000 $\Omega$ (nominal)	
Rate accuracy	±1% nominal	
ECG amplitudes		
Reference lead	Selectable, Lead II (default) or Lead I	
Settings	0.05 mV to 0.45 mV by 0.05 mV steps and 0.5 mV to 5 mV by 0.5 mV steps	
Accuracy (all performance waves and normal sinus R waves)		
Lead II	± 2 %	
All other leads	± 5 %	
Defib paddles	± 5 %	
Amplitude of ECG signals relative to amplitude setting (in percent)		
Lead II reference		
Performance waves and R wave detection		
Lead #	Ref. amp.	
1	70 %	
II	100 %	
	30 %	



V1	100 %	
V2	100 %	
V3	100 %	
V4	100 %	
V5	100 %	
V6	100 %	
Normal sinus waves:		
Lead #	Ref. amp.	
1	70 %	
11	100 %	
111	30 %	
V1	24 %	
V2	48 %	
V3	100 %	
V4	120 %	
V5	112 %	
V6	0 %	
Lead I reference		
Performance waves and R wave de	tection	
Lead #	Ref. amp.	
1	100 %	
II	150 %	
111	50 %	
V1	100 %	
V2	100 %	
V3	100 %	
V4	100 %	
V5	100 %	
V6	100 %	
Normal sinus waves:		
Lead #	Ref. amp.	
1	100 %	
11	150 %	
111	50 %	
V1	24 %	
V2	48 %	
V3	100 %	
V4	120 %	
V5	112 %	
V6	80 %	



ECG normal sinus		
Rates	10 BPM to 360 BPM in 1 BPM steps	
ECG high level output (BNC jack)		
Amplitude		
Range	0.5 V per mV of reference lead setting	
Accuracy	± 5 %	
Output impedance	50 Ω	
ECG on defibrillator input load		
Same as the Lead II amplitude but	limited to ± 4 mV	
ECG performance waves		
Square wave	2 Hz and 0.125 Hz	
Triangular wave	2 Hz and 2.5 Hz	
Sine waves	0.05, 0.5, 5, 10, 40, 50, 60, 100, 150, and 200 Hz	
Pulse	30 BPM and 60 BPM, 60 ms pulse width	
R-wave detection		
Waveform	Haver-triangle	
Amplitude	0.05 mV to 0.45 mV in 0.05 mV steps and 0.5 mV to 5 mV in 0.5 mV steps	
Rate	30, 60, 80, 120, 200, and 250 BPM	
Widths	8, 10, 12 ms, and 20 ms to 200 ms in 10 ms steps	
Accuracy	± (1 % setting + 1 ms)	





Noise immunity		
Wave	Sine	
Line frequency	50 Hz or 60 Hz (± 0.5 Hz)	
Amplitude		
Range	0.0 mV to 10 mV in 0.5 mV steps	
Accuracy	± 5 %	
Transvenous pacer pulse simulatio	n	
Widths		
Range	0.1 ms, 0.2 ms, 0.5 ms, 1 ms, and 2 ms	
Accuracy	± 5 % of setting	
Amplitudes		
Range	0 (off) and $\pm$ 2 mV, $\pm$ 4 mV, $\pm$ 6 mV, $\pm$ 8 mV, $\pm$ 10 mV, $\pm$ 12 mV, $\pm$ 14 mV, $\pm$ 16 mV, $\pm$ 18 mV, $\pm$ 20 mV, $\pm$ 50 mV, $\pm$ 100 mV, $\pm$ 200, $\pm$ 500, and $\pm$ 700 mV	
Accuracy	± (10 % setting + 0.2 mV)	
Amplitude of transvenous pacer pu	Ilse simulation signals relative to amplitude setting (in percent)	
Lead II reference		
Lead #	Ref. amp.	
1	67 %	
11	100 %	
111	33 %	
V1	67 %	
V2	67 %	
V3	67 %	
V4	67 %	
V5	67 %	
V6	67 %	
Lead I reference		
Lead #	Ref. amp.	
1	100 %	
II	150 %	
Ш	50 %	
V1	100 %	
V2	100 %	
V3	100 %	
V4	100 %	
V5	100 %	
V6	100 %	



Arrhythmia selections		
Pacer interactive (7000DP only)	Demand: 30 BPM to 360 BPM in 1 BPM steps	
	Asynchronous	
	Non-capture	
	Non-function	
	Threshold (interactive pacing simulation only): 10 mA to 250 mA in 10 mA steps	
Supraventricular	Artrial fibrillation course	
	Atrial fibrillation fine	
	Atrial flutter	
	Sinus arrhythmia	
	Missed beat	
	Atrial tachycardia	
	Paroxysmal atrial tachycardia (PAT)	
	Nodal rhythm	
	Supraventricular tachycardia Premature	
Premature	Atrial PAC	
	Nodal PNC	
	PVC1 left ventricle	
	PVC1 LV early	
	PVC1 LV R on T	
	PVC2 right ventricle	
	PVC2 RV early	
	PVC2 RV R on T	
	Multifocal PVCs	
Ventricular	PVCs 6/min	
	PVCs 12/min	
	PVCs 24/min	
	Freq multifocal	
	Trigeminy	
	Bigeminy	
	Pair PVCs	
	Run 5 PVCs	
	Run 11 PVCs	
	Monomorphic ventricular tachycardia: 120 BPM to 300 BPM in 5 BPM steps	
	Polymorphic ventricular tachycardia: 1 to 5	
	Ventricular fibrillation: coarse and fine	



### Conduction

	1° Block
	2° Block Type I
	2° Block Type II
	3° Block
	Right bundle branch block RBBB
	Left bundle branch block LBBB Transvenous Paced with select- able pacer spike amplitudes and widths
	Atrial 80 BPM
	Async 75 BPM
	Demand with frequent sinus beats
	Demand with occasional sinus beats
	AV sequential
	Non-capture
	Non-function

## Selections for all waves in group

Atrial pacer pulse		
Width	0.1, 0.2, 0.5, 1, 2 ms	
Polarity	+ or -	
Amplitude	0 (off), 2 to 20 (by 2), 50, 100, 200, 500, 700 mV	
Ventricular pacer pulse		
Width	0.1, 0.2, 0.5, 1, 2 ms	
Polarity	+ or -	
Amplitude	0 (off), 2 to 20 (by 2), 50, 100, 200, 500, 700 mV	
R-wave detection		
Rate	30, 60, 80, 120, 200, 250 BPM	
Width	8, 10, 12, 20 to 200 (by 10) ms	
Amplitude	0.05 to 0.45 (by 0.05), 0.5 to 5 (by 0.5) mV	

## Transcutaneous pacemaker analyzer technical specifications (7000DP only)

Test load Selections	
Defibrillator input	
Fixed load	50 Ω
Accuracy	$\pm$ 1 %, non-inductive (<2 $\mu\text{H})$
Power rating	10 defib pulses of 360 J every 5 minutes
Pacemaker input	
Variable load	50 $\Omega$ to 1500 $\Omega$ in 50 $\Omega$ steps
Accuracy	± 2 %, non-inductive (< 2 μH)
Power rating	5 W (average), 40 W (peak) @ 1000 Ω



Measurements	
Manufacturer specific algorithms	GE Responder (1500 and 1700)
	MDE 300 (Medical Data Electronics)
	Medtronic ERS/Physio Control LIFEPAK
	MRL (Medical Research Laboratory/Welch Allyn)
	Philips/Agilent/HP
	Schiller Medical
	ZOLL Medical (plus a general purpose Default Algorithm selection)
Current	
Range	4 mA to 250 mA
Accuracy	±1% of reading + 0.02 mA
Pulse rate	
Range	5 PPM to 800 PPM
Accuracy	± 0.5 % of reading + 0.1 PPM
Pulse width	
Range	1 ms to 100 ms
Accuracy	± 0.5 % of reading + 0.01 ms
Energy	
Range	1 μJ to 2 J
Accuracy	$\pm4$ % of reading + 10 $\mu J$
Demand and asynchronous mode	test
Input pacer pulse rates	30 PPM to 200 PPM
Energy	
Rate	10 BPM to 300 BPM in 1 BPM steps
Amplitude	1 mV
Underdrive rate	10 BPM minimum
Overdrive rate	300 BPM maximum
Sensitivity test	
Automatic interactive threshold de	tection
Compatible pacer rates	30 PPM to 120 PPM
ECG R wave	
Waveforms	Square, triangle, sine
Width	1 ms to 19 ms (by 1 ms), 20 ms to 95 ms (by 5 ms), 100 ms to 300 ms (by 25 ms)
Accuracy	±5% of setting
Amplitude	0.05 mV to 0.95 mV (by 0.05 mV), 1 mV to 5 mV (by 0.5 mV)
Accuracy	$\pm$ 5 % of setting



Refractory period tests		
Paced refractory period	20 ms to 500 ms	
Sensed refractory period	15 ms to 500 ms	
Accuracy	±1ms	
Pacer pulse rate	20 PPM to 200 PPM	
ECG		
Waveform	Triangle wave	
Pulse width	40 ms	
Amplitude	1 mV	





## Impulse 7010 Defibrillator Selectable Load Accessory

General specifications	
Maximum voltage	5000 V
Maximum continuous power	12 W, equivalent to 10 defib pulses of 360 J every 5 minutes
Inductance	< 2 μH, @25 Ω
	< 3 μΗ, @50 Ω
	< 4 $\mu H,$ @75 $\Omega$ and 100 $\Omega$
	< 5 μΗ, @125 Ω
	< 6 μΗ, @150 Ω
	< 7 μH, @175 Ω
	< 8 μH, @200 Ω
Operating Temperature	10 °C to 40 °C (50 °F to 104 °F)
Storage Temperature	-20 °C to 60 °C (-4 °F to 140 °F)
Humidity	10 % to 90 % non-condensing
Dimensions (WxDxH)	154 mm x 272 mm x 138.7 mm (6.07 in x 10.71 in x 5.46 in)
Weight (net)	1.54 kg (3 lb 6.2 oz)
Safety class	Complies with EN61010-1 2nd Edition, Class II product
Safety and EMC marks	
Warranty	Two-year extended warranty (no-cost extended warranty available after first-year calibration at any Fluke Biomedical autho- rized service center)
Calibration interval	One-year
Electrical specifications (for Load A	Accessory and Analyzer together)
Load settings	25 $\Omega,$ 50 $\Omega,$ 75 $\Omega,$ 100 $\Omega,$ 125 $\Omega,$ 150 $\Omega,$ 175 $\Omega,$ and 200 $\Omega\pm$ 1 $\%$
Accuracy	
Energy (all except pulsed biphasic)	2 % of reading + 0.1 J with 25, 75 $\Omega$ though 200 $\Omega$ loads, 1 % of reading + 0.1 J with 50 $\Omega$ load
Energy (pulsed biphasic)	2.5 % of reading + 0.3 J with 25, 75 $\Omega$ though 200 $\Omega$ loads, 1.5 % of reading + 0.3 J with 50 $\Omega$ load
Voltage	1 % of reading + 2 V with 25 $\Omega$ and 50 $\Omega$ loads, 2 % of reading + 2 V with 75 $\Omega$ through 200 $\Omega$ loads
Current	2 % of reading + 0.1 A with 25 $\Omega$ load, 1 % of reading + 0.1 A with 50 $\Omega$ through 200 $\Omega$ loads



## Ordering information

## Models/descriptions

2811928	Impulse 6000D Defibrillator Analyzer 120 V (US)
3077031	Impulse 6000D Defibrillator Analyzer (Schuko)
3077046	Impulse 6000D Defibrillator Analyzer (UK)
3077054	Impulse 6000D Defibrillator Analyzer (Japan)
3085270	Impulse 6000D Defibrillator Analyzer (Australia)
3085281	Impulse 6000D Defibrillator Analyzer (India)
2811919	Impulse 7000DP Defibrillator/ Transcutaneous Pacemaker Analyzer 120 V (US)
3077005	Impulse 7000DP Defibrillator/ Transcutaneous Pacemaker Analyzer 120 V (Schuko)
3077010	Impulse 7000DP Defibrillator/ Transcutaneous Pacemaker Analyzer 120 V (UK)
3077022	Impulse 7000DP Defibrillator/ Transcutaneous Pacemaker Analyzer 120 V (Japan)
3085296	Impulse 7000DP Defibrillator/ Transcutaneous Pacemaker Analyzer 120 V (Australia)
3085308	Impulse 7000DP Defibrillator/ Transcutaneous Pacemaker Analyzer 120 V (India)
3326874	TA-IMP7KDP Impulse 7000DP Defibrillator/Transcutaneous Pacemaker Ana- lyzer with test automation 120 V (US)
3326888	TA-IMP7KDP-01 Impulse 7000DP Defibrillator/Transcutaneous Pacemaker Analyzer with test automation 120 V (Schuko)
3326895	TA-IMP7KDP-02 Impulse 7000DP Defibrillator/Transcutaneous Pacemaker Analyzer with test automation 120 V (UK)
3326901	TA-IMP7KDP-03 Impulse 7000DP Defibrillator/Transcutaneous Pacemaker Analyzer with test automation 120 V (Japan)
3326912	TA-IMP7KDP-04 Impulse 7000DP Defibrillator/Transcutaneous Pacemaker Analyzer with test automation 120 V (Australia)
3326920	TA-IMP7KDP-05 Impulse 7000DP Defibrillator/Transcutaneous Pacemaker Analyzer with test automation 120 V (India)

## **Standard accessories**

1626219	USB Computer Communication Cable
3028662	Getting-Started Guide
Battery Eliminator	(country specific)
2814980	Carrying Case
3156262	Defib Paddle Contact Plates



## **Optional accessories**

3091370	Ansur Impulse 6000D/7000DP Plug-In
3065489	MedtronicERS/Physio-Control (FAST PATCH) (set of two): 4 mm defibrillator adapters
3065450	Kimberly Clark/R2 Darox MRL/MDE/ NK: 4 mm defibrillator adapters
3065438	Internal discharge paddle contacts (set of two)
3065477	Medtronic ERS/Physio-Control (QUIK PACE) (set of two): 4 mm pacer adapters
3065527	Zoll Medical NTP/PD1400: 4 mm pacer adapters
3065461	Medtronic ERS/Physio-Control (QUIK COMBO): 4 mm defib/pacer adapters
3065492	Philips/Agilent/HP (CODEMASTER Series-Round): 4 mm defib/pacer adapters
3065509	Philips/Agilent HEARTSTART FR2/ MRX: 4 mm defib/pacer adapters
3065511	Zoll PD-2200 Multi-Function PDSeries, M-Series, M-Series CCT, AED PRO and AED Plus <sup>™</sup> defib/pacer adapters
3065423	GE Marquette (RESPONDER1500/1700 Series) (set of two): 4 mm defib/ pacer adapters
3158544	Impulse 7010 Defibrillator Selectable Load Accessory





### **About Fluke Biomedical**

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology qualityassurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-0 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

### Fluke Biomedical regulatory commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 and ISO 13485 medical device certified and our products are:

- CE Certified, where required
- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required
- NRC Compliant, where required

### Fluke Biomedical.

Trusted for the measurements that matter.

#### Fluke Biomedical

28775 Aurora Road Cleveland, OH 44139 U.S.A.

For more information, contact us at: (800) 850-4608 or Fax (440) 349-2307 Email: sales@flukebiomedical.com Web access: www.flukebiomedical.com

©2013, 2019 Fluke Biomedical. Specifications subject to change without notice. Printed in U.S.A. 10/2019 3393310d-en

Modification of this document is not permitted without written permission from Fluke Corporation.

FLUKE ®