




INSTRUCTIONS FOR USE  
BEDIENUNGSANLEITUNG  
INSTRUCTIONS D'UTILISATION  
INSTRUCCIONES DE USO  
ISTRUZIONI PER L'USO  
GEBRUIKSAANWIJZING  
BRUKSANVISNING  
INSTRUÇÕES DE USO  
使用说明

# RaySafe i3

ENGLISH.....	3
DEUTSCH .....	4
FRANÇAIS.....	5
ESPAÑOL.....	6
ITALIANO .....	7
NEDERLANDS.....	8
SVENSKA.....	9
PORTUGUÊS.....	10
简体中文.....	11
APPENDIX.....	12
EU DECLARATION OF CONFORMITY.....	13

© 2019.05 Unfors RaySafe 5001106-2.10

All rights are reserved. Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without the prior written consent of the copyright owner.

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## Intended use

RaySafe i3 is designed for personal use on hospital staff. It measures scattered radiation in examination rooms using a dosimeter, and shows a graphical visualization of dose rate in real time on a display. (Figure “1.” on page 2.)

## Usage

The RaySafe i3 system needs to be installed and initiated before use. This can be done by a service technician, using the instructions in the printed *Installation and Service Manual*, delivered with the system.

Use the dosimeter on the outside of your X-ray protection garment. Front forward, clip towards your body.

The Real-time Dosimeter (A in Figure “2.” on page 2) measures scattered radiation and stores measurements. When a dosimeter is exposed to radiation while in range of a Real-time Display (B in Figure “2.” on page 2), measurements are transferred from dosimeter to display via radio. The display shows dose rate bars in real time in green, yellow or red, from low to high dose rate. You can tap on the rows on the display to see more details and dose history.

To view dose history stored in the dosimeter, connect the dosimeter to a computer running Dose Viewer (C in Figure “2.” on page 2). See the *Dose Viewer manual*, found in the help menu in Dose Viewer.

If you require further information about product training, contact your local RaySafe representative.

**CAUTION!** *The Real-time Dosimeter measures scattered radiation and should neither be positioned in the primary X-ray field nor be used for diagnostic measurements.*

**WARNING!** *RaySafe i3 does not replace any X-ray protection equipment.*

**WARNING!** *RaySafe i3 is not sterile.*

**WARNING!** *RaySafe i3 is not a replacement for legal dosimeters.*

## Maintenance

The Real-time Dosimeter has a replaceable battery. When a red light is blinking on the side of the dosimeter, the battery needs to be replaced. A service technician can do this, using the instructions in the service manual.

The Real-time Display has no serviceable parts.

For cleaning of dosimeter and display, use a damp cloth and mild detergent.

**WARNING!** *Do not modify the Real-time Display.*

## Specifications

### Physical

	Real-time Dosimeter	Real-time Display
Weight	34 g (1.2 oz)	1.2 kg (42 oz)
Dimensions	40 × 58 × 17 mm (1.6 × 2.3 × 0.7 in)	300 × 250 × 60 mm (11.8 × 9.8 × 2.4 in)
Power	Replaceable battery: CR2450 Li/MnO <sub>2</sub>	External supply: Astec DPS53-M, GTM96600-6012-T3 or FW7556M/12

### Radiology

Operational quantity	H <sub>p</sub> (10)
Energy dependence	< 25 % (N-series, 40 – 150 kV)
Temperature dependence	< 5 % (18 – 26 °C) < 25 % (15 – 18 °C, 26 – 35 °C)
Response time	< 1 s (> 100 μSv/h), < 5 s (< 100 μSv/h)
Integration interval	1 s
Detection limit	< 30 μSv/h
Dose rate uncertainty (continuous radiation)	10 % or 10 μSv/h (40 μSv/h – 150 mSv/h), 20 % (150 – 300 mSv/h), 40 % (300 – 500 mSv/h)
The measured dose rate is monotonically increasing up to 1 Sv/h.	
Dose reproducibility	10 % or 1 μSv
Maximum lifetime dose	10 Sv
Typical angular response (relative to 0°)	Vertical: <i>Figure “3.” on page 2</i> Horizontal: <i>Figure “4.” on page 2</i>
Typical energy response	<i>Figure “5.” on page 2</i>
Real-time Display output unit	Sv or rem

## Radio

	Real-time Dosimeter	Real-time Display
Frequency	868.3 MHz or 918.3 MHz, depending on purchased configuration (see label on product)	
Transmitter power output	3.2 mW	1 mW
Communication range	Typically 5 – 15 m, depending on the local environment	

## Climatic conditions

Temperature	+15 – +35 °C (for operation) –25 – +70 °C (for storage)
Temperature rate of change	< 0.5 °C/minute
Relative humidity	20 – 80 % (for operation) 5 – 95 % (for storage)
IP code	Real-time Dosimeter: IP20 Real-time Display: IP21

## For safe usage

Always power off the dosimeter before air transport, to make sure that the radio communication is completely turned off.

If there is any risk of liquid spilling on the Real-time Display, turn it off by disconnecting the power cord.

## Security and privacy

Dosimeter serial number as well as name and settings entered via Dose Viewer are sent unencrypted via radio to the Real-time Display.

## Symbols on labels

Symbol	Meaning
	Manufacturer
	Date of manufacture
	Article number
	Serial number
	Consult the instructions for use.
	Product complies with essential requirements applicable to European directives
	Electrical safety compliance marking, together with listing number (“Appendix” on page 12)
	Radio certification number for US market (“Appendix” on page 12)
IC	Radio certification number for Canadian market (“Appendix” on page 12)
	Radio certification symbol and number for Brazilian market (“Appendix” on page 12)
	In the European Union the WEEE directive and label indicates that when discarding this product, it must be sent to separate collection facilities for recovery and recycling. For advice and information, contact RaySafe.

## Learn more

- Real-time Dosimeter quick guide under the lid of the dosimeter packaging
- Dose Viewer manual via the Help button in Dose Viewer
- Installation and Service manual delivered with the RaySafe i3 system

## Contact information

For technical support, contact [technicalsupport@raysafe.com](mailto:technicalsupport@raysafe.com).

Instructions and contact information can also be found at <http://www.raysafe.com>

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## MET

MET Laboratories, Inc. Certification covers UL 60950-1/CSA C22.2 No. 60950-1. MET Laboratories, Inc. has not evaluated this device for reliability or efficacy of its intended functions.

## FCC

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

*NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.*

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This portable device with its antenna is excluded from SAR testing since the output power is below the SAR test exclusion threshold. The antenna used for this device must not be co-located or operating in conjunction with any other antenna or transmitter.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## IC/ISED

This Device complies with Industry Canada License-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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**This declaration of conformity is issued under the sole responsibility of the manufacturer.**

Product name: RaySafe i3 System 868,3 MHz  
Product Model Number: 160608x, where x is a digit with no impact on the declaration of conformity.

Product Options/Accessories: Parts as described in the accompanying documentation

**The object of the declaration described above is in conformity with:**

- Directive **2014/53/EU** of the European Parliament and of the Council of 16 April 2014 on the harmonisation of the laws of the Member States relating to the making available on the market of radio equipment and repealing Directive 1999/5/EC.
- Directive **2011/65/EU** of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

The product(s) listed above have been tested in a typical configuration as described in the Manufacturer's accompanying documentation, and are fully compliant with the standards listed below. Additionally the products listed above have been designed, manufactured, tested, and found to be compatible with the devices and accessories described by the manufacturer in the devices accompanying documentation.

Date of Issue: 2017-Sep-1

Signature:  
(signed for and on  
behalf of RaySafe)

Name:  
Title:

Jesper Bengtsson  
Manager RA/QA

Place of issue: Billdal, Sweden

The device as mentioned above to which this EU Declaration of conformity is applicable, has been developed, designed and found to be in conformity pursuant to the European harmonized standards specified hereafter

<b>Harmonized Standard:</b>	<b>Title:</b>
EN60950-1:2006 /A11:2009/A1:2010/ A12:2011/A2:2013	Information technology equipment - Safety – Part 1: General requirements
ETSI EN 301 489-1 V2.1.1	Electromagnetic Compatibility (EMC)= standard for radio equipment and services; Part 1: Common technical requirements; Harmonised Standard covering the essential requirement of article 3.1(b) of Directive 2014/53/EU and the essential requirements of article 6 of Directive 2014/30/EU
ETSI EN 301 489-3 V2.1.1	Electromagnetic Compatibility (EMC) standard for radio equipment and services; Part 3: Specific conditions for Short-Range Devices (SRD) operating on frequencies between 9 kHz and 246 GHz; Harmonised Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU
ETSI EN 300 220-1 V3.1.1	Short Range Devices (SRD) operating in the frequency range 25 MHz to 1 000 MHz; Part 1: Technical characteristics and methods of measurement
ETSI EN 300 220-2 V3.1.1	Short Range Devices (SRD) operating in the frequency range 25 MHz to 1 000 MHz; Part 2: Harmonised Standard covering the essential requirements of article 3.2 of Directive 2014/53/EU for non specific radio equipment
EN 50581:2012	Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

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Name / address manufacturer: Unfors RaySafe AB