

# Manual Supplement

Manual Title:	Impulse 6000D/7000DP Getting Started	Supplement Issue:	7
Part Number	3028662	Issue Date:	4/13
Print Date:	August 2007	Page Count:	9
Revision/Date:			

---

This supplement contains information necessary to ensure the accuracy of the above manual.

**FLUKE**®

**Biomedical**

## **Change #1, 43260, 43305, 45041, 44082, 46470, 46697, 47714, 48007, 49680**

Replace pages 15 through 19, with the following:

### ***Defibrillator Analyzer Specifications***

#### ***Energy Output Measurement***

**Compatible Defibrillator Waveshapes** ..... Lown, Edmark, Trapezoidal, DC Bi-phasic, and AC Pulsed Bi-phasic

*Note*

*AC Pulsed Biphasic waveform has not been approved in the United States.*

**Autoranged Measurement**..... 0.1 to 600 J

#### **Accuracy**

0.1 to 360 J .....  $\pm(1\% \text{ of reading} + 0.1 \text{ J})$

360 to 600 J .....  $\pm(1\% \text{ of reading} + 0.1 \text{ J})$ , typical

*Note*

*For Pulsed Bi-Phasic defibrillator, specified accuracy is  $\pm(1.5\% \text{ of reading} + 0.3 \text{ J})$  on both ranges.*

#### **Load resistance**

Resistance ..... 50  $\Omega$

Accuracy .....  $\pm 1\%$ , non-inductive (<2  $\mu\text{H}$ )

**Pulse trigger level** ..... 20 V

#### **Pulse width**

Range ..... 1.0 to 50.0 ms

Accuracy .....  $\pm 0.1$  ms

**Voltage**

Range ..... 20 to 5000 V

Accuracy .....  $\pm(1\%$  of reading + 2 V)

**Current**

Range ..... 0.4 to 100.0 A

Accuracy .....  $\pm(1\%$  of reading + 0.1 A)

**Tilt (biphasic and pulsed biphasic)**

Range ..... 1 % to 99 %

Accuracy .....  $\pm 1$  digit

**Interphase delay (biphasic and pulsed biphasic)**

Range ..... 0.1 ms to 9.9 ms

Accuracy .....  $\pm 0.1$  ms

**Frequency (pulsed biphasic only)**

Range ..... 2000 Hz to 8000 Hz

Accuracy .....  $\pm 1\%$  of reading

**Duty cycle (pulse biphasic only)**

Range ..... 1 % to 99 %

Accuracy .....  $\pm 1$  digit

**Sample rate** ..... 250 kHz (4  $\mu$ s sample)

**Maximum Average Power** ..... 12 W, equivalent to 10 defib pulses of 360 J every 5 minutes

**Scope Output**

Autorange ..... 2000:1, 400:1 and 80:1: dependant on the range

Waveform Playback

Output ..... BNC

Output impedance..... 50 Ω

Amplitude Accuracy ..... ±5 %

**Charge Time Measurement**

Range ..... 0.1 to 100.0 s

Accuracy ..... ±0.05 s, typical

**Synchronization Test (Elective Cardioversion)**

Delay Time Measurement

Timing window ..... ECG R-wave peak to the defib pulse peak

Range ..... -120 to +380 ms; measures timing from 120 ms prior to the R-wave peak to up to 380 ms following the R-wave peak.

Resolution ..... 1 ms

Accuracy ..... ±1 ms

ECG waves

Normal Sinus Rhythm (NSR)..... 10 to 180 (by 1) BPM

Atrial fibrillation..... Coarse and fine

Monomorphic Ventricular Tachycardia ..... 120 to 240 (by 5) BPM

Asystole ..... Flat line

**Automated Defibrillator Test ECG Waves**

Normal Sinus ..... 10 to 300 (by 1) BPM

- Ventricular Fibrillation ..... Coarse and fine
- Monomorphic Ventricular Tachycardia ..... 120 to 300 (by 5) BPM
- Polymorphic Ventricular Tachycardia ..... 5 types
- Asystole ..... Flat line

**ECG Waves**

**ECG General**

- Lead configuration ..... 12-lead simulation. RA, LL, LA, RL, V1-6 with independent outputs
- Lead to lead impedance ..... 1000 Ω (nominal)
- Rate accuracy ..... ±1 % of nominal

**ECG Amplitudes**

- Reference Lead ..... Selectable, Lead II (default) or Lead I
- Settings ..... 0.05 to 0.45 (by 0.05) mV  
0.5 to 5.0 (by 0.5) mV

Accuracy (All Performance waves and Normal Sinus R waves)

- Lead II ..... ±2 %
- All other leads ..... ±5 %
- Defib paddles ..... ±5 %

Amplitude of ECG signals relative to amplitude setting (in percent)

Lead II reference

Performance waves and R wave detection

Lead #	I	II	III	V1	V2	V3	V4	V5	V6
Ref. Amp.	70 %	100 %	30 %	100 %	100 %	100 %	100 %	100 %	100 %

Normal Sinus waves

Lead #	I	II	III	V1	V2	V3	V4	V5	V6
Ref. Amp.	70 %	100 %	30 %	24 %	48 %	100 %	120 %	112 %	80 %

Lead I reference

Performance waves and R wave detection

Lead #	I	II	III	V1	V2	V3	V4	V5	V6
Ref. Amp.	100 %	150 %	50 %	100 %	100 %	100 %	100 %	100 %	100 %

Normal Sinus waves

Lead #	I	II	III	V1	V2	V3	V4	V5	V6
Ref. Amp.	100 %	150 %	50 %	24 %	48 %	100 %	120 %	112 %	80 %

**ECG Normal Sinus**

Rates ..... 10 to 360 (by 1) BPM

**ECG High Level Output (BNC Jack)**

Amplitude

Range ..... 0.5 V per mV of reference lead setting

Accuracy ..... ±5 %

Output Impedance ..... 50 Ω

**ECG on Defibrillator Input Load**

Same as the LEAD II amplitude but limited to ±4 mV

**ECG Performance Waves**

Square wave ..... 2.0 and 0.125 Hz

Triangular wave ..... 2.0 and 2.5 Hz



Accuracy .....  $\pm(10\% \text{ of setting} + 0.2 \text{ mV})$

Amplitude of Transvenous Pacer Pulse Simulation signals relative to amplitude setting (in percent)

Lead II reference

Lead #	I	II	III	V1	V2	V3	V4	V5	V6
Ref. Amp.	67 %	100 %	33 %	67 %	67 %	67 %	67 %	67 %	67 %

Lead I reference

Lead #	I	II	III	V1	V2	V3	V4	V5	V6
Ref. Amp.	100 %	150 %	50 %	100 %	100 %	100 %	100 %	100 %	100 %

**Arrhythmia Selections**

Pacer Interactive (Transcutaneous pacer, Impulse 7000DP only)

Demand ..... 30 to 360 (by 1) BPM

Asynchronous

Non-Capture

Non-Function

Threshold (Interactive pacing simulation only) 10 to 250 (by 10) mA

Supraventricular

Atrial Fibrillation Coarse

Atrial Fibrillation fine

Atrial Flutter

Sinus Arrhythmia

Missed Beat

Atrial Tachycardia



## Change #2

On page 1, under ***Intended Use***, replace the current paragraph with following:

The Analyzer is used to determine if defibrillators and transcutaneous pacemakers are performing within their performance specifications through measurement of energy output.

## Change #3, 48038

On page 22, under **Pacemaker Input**, change:

From: Accuracy .....±1 % non-inductive (<2 μH)

To: Accuracy.....±2 % non-inductive (<2 μH)

## Change #4


On page 3 add the following Caution:

### Caution

**To avoid damage to the product or defibrillator under test, do not use adhesive pads on the Analyzer. For defibrillators with metal paddles, use Fluke Biomedical electrode adapters.**

## Change #5, 63282, 65427

On page 2, delete CAT I from the **Symbols** table and add:

	Conforms to relevant South Korean EMC Standards.
---	--

On page 14, under **General Specifications** remove **Safety Standards** and **Electromagnetic Compatibility Standards (EMC)** and add:

Safety .....	IEC 61010-1: CAT II 300 V, Pollution Degree 2
Altitude .....	2000 m
IP Rating.....	IEC 60529: IP20
Electromagnetic Environment.....	IEC 61326-1: Basic
Electromagnetic Compatibility .....	Applies to use in Korea only, Class A Equipment (Industrial Broadcasting & Communication Equipment) <sup>[1]</sup>

[1] This product meets requirements for industrial (Class A) electromagnetic wave equipment and the seller or user should take notice of it. This equipment is intended for use in business environments and is not to be used in homes.