

Biomedical

Masimo Rainbow SpO₂ technology and testing

Application Note



Introduction


Masimo Rainbow is the first technology to non-invasively measure blood constituents and fluid responsiveness that previously required invasive procedures. Besides measuring arterial oxygensaturation level (SpO_2) and pulse rate, Rainbow technology measures following parameters as well: Hemoglobin (SpHb), Carboxyhemoglobin (SpCO) and Methemoglobin (SpMet).

Traditional pulse oximeters use wavelengths of two lights (infrared light at 940 nm and red light at 660 nm), which are absorbed differently by either oxyhemoglobin or the reduced hemoglobin to measure SpO₂ concentration. Two-wavelength oximeters cannot measure total hemoglobin or dyshemoglobins. Rainbow sensor technology uses more than 7 wavelengths of light to acquire blood constituent data based on light absorption

When used with the Masimo Rainbow Test Sensor, ProSim Vital Sign Simulator is the first and only simulator capable of verifying the performance of Masimo Rainbow Oximeters.

How to test Masimo Rainbow

Attach the Rainbow Test sensor to the artificial finger of ProSim 8 as shown in the figure below. Place the sensor with the LEDs on the bottom of the artificial finger. While placing the sensor on the artificial finger, monitor the signal indicator along the bottom of the product display. Adjust the sensor on the finger for maximum signal strength.





Biomedical

2. Select Masimo Rainbow from Type selection:



3. The SpO₂ screen shows three more parameters than what is shown for other types of sensors: SpMet, SpCO, and SpHb.



- a. SpMet, SpCO, and SpHb cannot be set through the ProSim. The special Masimo Rainbow cable sets them based on the measured SpO₂ percent.At 100 %, SpMet = 0 %, SpCO = 0%, and SpHb = 25 g/dl.
 - i. SpCO: -1 % change in SpO₂ = +0.7 % change in SpCo
 - ii. SpMet: -1 % change in $SpO_2 = +0.3$ %, SpMet
 - iii. SpHb: -1 % change in $SpO_2 = -0.5$ g/dl change in SpHb
 - iv. SpHb does not change for values of SpO₂ above 90 %.

About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-0 accredited laboratory, Fluke Biomedical also offers the best

in quality and customer service for all your equipment calibration needs. Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment
As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 and ISO 13485 medical device certified and our products are:

CE Certified, where required

NIST Traceable and Calibrated

- UL, CSA, ETL Certified, where requiredNRC Compliant, where required

Fluke Biomedical.

Better products. More choices. One company.

Fluke Biomedical

6045 Cochran Road Cleveland, OH 44139-3303 U.S.A.

Fluke Biomedical Europe

Science Park Eindhoven 5110 5692EC Son, The Netherlands

For more information, contact us:

In the U.S.A. (800) 850-4608 or Fax (440) 349-2307 In Europe/M-East/Africa +31 40 267 5435 or Fax +31 40 267 5436 From other countries +1 (440) 248-9300 or Fax +1 (440) 349-2307 Email: sales@flukebiomedical.com Web access: www.flukebiomedical.com

©2012 Fluke Biomedical. Specifications subject to change without notice. Printed in U.S.A. 5/2012 4207999A_EN

Modification of this document is not permitted without written permission from Fluke Corporation.