

# Manual Supplement

Manual Title: SPOT Light Users  
Part Number: 4151274  
Print Date: February 2012  
Revision/Date: 1

Supplement Issue: **3**  
Issue Date: 8/17  
Page Count: 1

---

---

This supplement contains information necessary to ensure the accuracy of the above manual.


### Change #1, 63286, 218

On pages 14, under **Specifications** remove **Battery Charger, Safety Standards, Certifications** and **Electromagnetic Compatability** and replace with:

- Safety** ..... IEC 61010-1: 300 V CAT II, Pollution Degree 2
- IP Rating** ..... IEC 60529: IP20
- Electromagnetic Environment** ..... IEC 61326-1: Basic
- Battery Charger** ..... 100 V to 240 V, 50/60 Hz input, 6V/3 A output. For best performance, the battery charger must be connected to a properly grounded ac receptacle.

### Change #2, 456, 458

On page 9, add the following to the end of the page:

By default, the unit automatically turns off if inactive for 5 minutes. To change this setting push and hold  for 3 seconds when the 'Test' field is highlighted. This shows an auto-off option. Select Y for Yes or N for No.

To restore the display to English, turn off SPOTLight, then push and hold  and turn on SPOTLight.

### Change #3, 475

On page 5, replace the AC/DC Power Supply with:

AC/DC Power Supply	4766589
--------------------	---------