

Manual Supplement

Manual Title:	VT MOBILE Getting Started	Supplement Issue:	4
Part Number:	2544892	Issue Date:	6/16
Print Date:	January 2006	Page Count:	3
Revision/Date:			

This supplement contains information necessary to ensure the accuracy of the above manual.

Change #1

On page 28, under *High-Flow Port*, change:

From: **Maximum Flow Rate**
(absolute value): 200 lpm

To: **Maximum Flow Rate**
(absolute value): 220 lpm

From: **Operating Flow Range:** ± 150 lpm

To: **Operating Flow Range:** ± 200 lpm

Change #2

On page's 2 and 3, Tables 2 and 3, replace Oxygen Sensor with:

Oxygen Sensor	2248801
---------------	---------

Change #3, 63288, 64415

On page vii, under the **Warnings** add:

- **Do not use the Product around explosive gas, vapor, or in damp or wet environments.**
- **Use the Product only as specified, or the protection supplied by the Product can be compromised.**

On page 26, add the following under **General Specifications**:




- Safety IEC 61010-1: CAT II 300 V, Pollution Degree 2
- Operating Altitude..... 2000 m
- Electromagnetic Environment IEC 61326-1: Basic
- Electromagnetic Compatibility Applies to use in Korea only. Class A Equipment (Industrial Broadcasting & Communication Equipment) ^[1]

[1] This product meets requirements for industrial (Class A) electromagnetic wave equipment and the seller or user should take notice of it. This equipment is intended for use in business environments and is not to be used in homes.

On page 31, under *Peak Expiratory Flow (PRF)* replace Accuracy with:

Accuracy: ± 3 % or 2 lpm

On page 33, add the following to the **Symbols** table:

	Conforms to relevant North American Safety Standards.
	Conforms to relevant South Korean EMC Standards.
	Conforms to relevant Australian Standards

Change #4, 207, 304

On page 3, Table 3, replace the Power Adapter, with:

Power Adapter, Universal (USA and International)	4760105
--	---------

On page 26, under ***External Power Supply***, replace the **Output Current** with:

Output Current: 2 A