

VT PLUS HF

Gas-Flow Analyzer

Technical Data



The VT PLUS HF is Fluke Biomedical's premier general-purpose gas-flow analyzer. In addition, special display modes and bi-directional flow make it perfect for fully and efficiently testing both conventional mechanical ventilators and high-frequency ventilators. EC.6.20 now requires 100 % completion of scheduled life-support device preventive maintenance every year, and VT PLUS HF can help meet those requirements. Multiple special-function tests make troubleshooting quick and efficient.

VT PLUS HF has the capability to measure either high- or low-flow and pressure, replacing the need for gauges and flow meters. It measures 21 ventilator parameters and can display all of them on one screen. Results can be printed directly from the unit or from a PC with included Windows-compatible software. VT PLUS HF also has onboard graphing capability and shows the minimum, maximum, average, and absolute measurement for all parameters.

Learning to use the VT PLUS HF is simple. Technicians control the unit using the VT PLUS HF user-friendly command system, or, if they're familiar with the RT-200, they can switch to a special control mode that uses RT-200-style commands.

VT PLUS HF can be operated with a variety of precision test lungs to ensure that ventilators are tested to manufacturers' specifications and clinical expectations with a fully NIST-traceable testing system.

Key features

- Bi-directional flow, pressure, volume, and oxygen concentration, and pressure measurements
- Low- and high-pressure, and flow measurement capability
- Special HF mode—up to 900 BPM (15 Hz)
- RS-232 and printer ports
- Included Windows-compatible graphics software
- All 21 ventilator parameters displayed at once on one screen
- Operation by user-friendly VT PLUS HF command mode or special RT-200 command mode
- Minimum, maximum, average, absolute, and graph for all parameters
- Multiple special-function tests for efficient troubleshooting

Optional features

- Operation with a variety of precision test lungs available from Fluke Biomedical to complete a fully NIST-traceable ventilator testing system

Specifications

| | | |
|-----------------------------|---|---|
| Power | 100 V ac to 240 V ac, 50/60 Hz | |
| Maximum over-voltage | 264 V ac | |
| Power consumption | < 132 V A | |
| Fuse rating | 0.5 A, slow blow | |
| Display | 320 x 240 LCD with CFL backlight | |
| Viewing area | 10.1 cm x 8.2 cm (3 in x 4 in), blue on white background | |
| Operational modes | Manual mode for simple tests or troubleshooting; computer-control mode, using RS-232 serial port for special applications; use of VT PLUS HF with VT for Windows software for recording graphs and logging data to a computer | |
| Output ports | RS-232 serial port, and parallel-printer port | |
| Oxygen measurement | | |
| Range | 0 % to 100 % | |
| Accuracy | ± 2 % FSO | |
| Resolution | 0.1 % O ₂ | |
| Transducer location | Internal | |
| Gas | | |
| Compatibility | Air, O ₂ , CO ₂ , N ₂ , N ₂ O, He, mixtures, or user-defined | |
| Reference units | ATP, STPDO, STPD21 and BTPS | |
| Test parameters | | |
| Continuous flow | Low flow | ± (2 % of reading and 1 % of range) |
| | High flow | ± (2 % of reading and 1 % of range) |
| Volumetric flow | Low-flow | |
| | Flow range | -25 lpm to 25 lpm |
| | Accuracy | ± 2 % of reading or ± 1 % of range, whichever is greater |
| | Frequency response | > 25 Hz or t ₁₀₋₉₀ < 40 ms, whichever is greater |
| | Low-flow dropout | 0.01 lpm |
| | Breath-detect threshold | 0.5 lpm |
| | Maximum-flow rate | 50 lpm |
| | Volume range | > ± 60 l |
| | Sample rate | 100 Hz |
| | Resolution | 0.01 lpm flow > 1 lpm; 0.001 lpm flow < 1 lpm |
| | Dynamic resistance | < 2.5 cmH ₂ O @ 5 lpm |
| | Fittings | 15 mm OD, 1:40 conical male; 0.25 in NPT ID per ASTM F-1054 aluminum with black anodized finish |
| | Notes: | |
| | <ul style="list-style-type: none"> • Tidal-volume accuracy: ± 3 % of reading or ± 2 ml, whichever is greater • Volume accuracy tested to 1 liter • Flow accuracy is specified for dry air or oxygen • Below 3.0 lpm, measurement accuracy is obtained by allowing the VT PLUS HF to fully warm up or manually zeroing before reading or documenting measurement | |
| | High-flow | |
| Flow range | -300 lpm to 300 lpm | |
| Accuracy | ± 2 % of reading or ± 2 % of range, whichever is greater | |

| | | | |
|--|---|--|--|
| Volumetric flow cont. | Frequency response | > 25 Hz | |
| | High-flow dropout | 25 lpm | |
| | Breath-detect threshold | 2 lpm | |
| | Maximum-flow rate | 500 lpm | |
| | Volume range | > ± 60 l | |
| | Dynamic resistance | < 2 cmH ₂ O @ 60 lpm | |
| | Sample rate | 100 Hz | |
| | Resolution | 0.01 lpm | |
| | Fittings | 22 mm OD, 1:40 conical male; 15 mm ID, 1:40 conical female per ASTM F-1054 aluminum with black anodized finish | |
| | Notes: | | |
| | • Tidal-volume accuracy: ± 3 % of reading or ± 10 ml, whichever is greater | | |
| | • Volume accuracy tested to 7 liters | | |
| | • Flow accuracy is specified for dry air or oxygen | | |
| | Low-pressure | | |
| | Range | ± 500 mmHg (10 psi) | |
| | Accuracy | ± 0.8 % of reading or ± 1.5 mmHg, whichever is greater | |
| | Frequency response | > 10 Hz | |
| | Resolution | 0.1 mmHg | |
| | Fittings | Luer lock, stainless steel | |
| | Maximum applied pressure | 60 psi | |
| | Sample rate | 100 Hz | |
| | Operating pressure | 30 psi | |
| | Note: Fluid pressure may be applied to the positive port; however, fluids should be kept from entering the pressure port by using a suitable length of connection tubing | | |
| | High-pressure | | |
| | Maximum applied pressure | 150 psi | |
| | Range | ± 100 psi | |
| | Accuracy | ± 1 % of reading or ± 0.3 psig, whichever is greater | |
| | Frequency response | > 10 Hz | |
| | Resolution | 0.1 psi | |
| | Sample rate | 100 Hz | |
| | Fittings | DISS connector, stainless steel | |
| | Airway-pressure | | |
| Maximum applied pressure | 20 psi | | |
| Range | ± 120 cmH ₂ O | | |
| Accuracy | ± 0.75 % of reading or ± 0.5 cmH ₂ O, whichever is greater | | |
| Frequency response | > 25 Hz or t10-90 < 40 ms, whichever is greater | | |
| Resolution | 0.1 cmH ₂ O | | |
| Sample rate | 100 Hz | | |
| Fittings | Internally connected at the transducer distal end | | |
| Note: Airway pressure is internally tapped off the proximal-flow sensor port, which is the port closest to the exhaust port on the VT PLUS HF | | | |

| Ventilator parameter | | |
|---|------------------------|--|
| Inspiratory and expiratory tidal volume | Resolution | 0.1 ml |
| | Range | As specified in high-flow/low-flow specification |
| | Accuracy | As specified in high-flow/low-flow specification |
| Expiratory minute volume | Resolution | 0.001 lpm |
| | Range | 0 L to 60 L |
| | Accuracy | ± 3 % |
| Breath rate | Resolution | 0.1 BPM |
| | Range | 0.5 BPM to 150 BPM |
| | Accuracy | ± 1 % |
| Inspiratory-to-expiratory time ratio (I:E ratio) | Resolution | 0.01 |
| | Range | 1:200 to 200:1 |
| | Accuracy | ± 2 % or ± 0.1 s |
| Inspiratory time | Resolution | 0.01 s |
| | Range | 0 s to 60 s |
| | Accuracy | ± 1 % or ± 0.02 s |
| Expiratory time | Resolution | 0.01 s |
| | Range | 0 s to 90 s |
| | Accuracy | ± 1 % or ± 0.01 s |
| Peak inspiratory pressure | Resolution | 0.1 cmH ₂ O |
| | Range | ± 120 cmH ₂ O |
| | Accuracy | ± 3 % or ± 1 cmH ₂ O |
| Inspiratory pause pressure | Resolution | 0.1 cmH ₂ O |
| | Range | ± 120 cmH ₂ O |
| | Accuracy | ± 3 % or ± 1 cmH ₂ O |
| Mean airway pressure | Resolution | 0.1 cmH ₂ O |
| | Range | ± 80 cmH ₂ O |
| | Accuracy | ± 3 % or ± 0.5 cmH ₂ O |
| Positive end-expiratory pressure (PEEP) | Resolution | 0.1 cmH ₂ O |
| | Range | -5 cmH ₂ O to 40 cmH ₂ O |
| | Accuracy | ± 3 % or ± 0.5 cmH ₂ O |
| Lung compliance | Resolution | 0.1 ml/cmH ₂ O |
| | Range | 0 ml/cmH ₂ O to 150 ml/cmH ₂ O |
| | Accuracy | ± 5 % or ± 5 ml/cmH ₂ O |
| | Inspiratory pause time | > 0.5 s |
| Inspiratory hold time | Resolution | 0.01 s |
| | Range | 0 s to 60 s |
| | Accuracy | ± 1 % or ± 0.1 s |
| Expiratory hold time | Resolution | 0.01 s |
| | Range | 0 s to 90 s |
| | Accuracy | ± 1 % or ± 0.1 s |

| | | |
|---|--|--------------------|
| Peak expiratory flow | Resolution | 0.01 lpm |
| | Range | 0 lpm to 300 lpm |
| | Accuracy | ± 3 % or ± 2 lpm |
| Peak inspiratory flow | Resolution | 0.01 lpm |
| | Range | 0 lpm to 300 lpm |
| | Accuracy | ± 3 % or ± 2 lpm |
| Flow bias | Resolution | 0.01 lpm |
| | Range | 0 lpm to 30 lpm |
| | Accuracy | ± 2 % or ± 0.5 lpm |
| | Expiratory pause time | > 0.5 s |
| Operating environment conditions | | |
| Temperature range | 10 °C to 40 °C | |
| Ambient humidity | 0 % to 80 % non-condensing to 31 °C, decreasing to 50 % at 40 °C | |
| Barometric pressure | 8 psig to 18 psig | |
| Storage environment conditions | | |
| Temperature range | -25 °C to -50 °C | |
| Humidity | 0 to 95 % non-condensing | |
| Dimensions (WxDxH) | 25.4 cm x 25.4 cm x 12.7 cm (10 in x 10 in x 5 in) | |
| Weight | 4.53 kg (10 lb) | |

Ordering information

Models

VT+HF-US120 United States, 120 V
VT+HF-AUS250V Australia, 250 V
VT+HF-SHK250V Schuko, 250 V
VT+HF-BRAZ250 Brazil, 250 V
VT+HF-UK250V United Kingdom, 250 V

Premium precision ventilator test kits

(VT PLUS HF Gas-Flow Analyzer; and ACCU LUNG portable precision test lung)

VT+HF/ACCULUNG-US United States
VT+HF/ACCULUNG-AUS Australia
VT+HF/ACCULUNG-SHK Schuko
VT+HF/ACCULUNG-BRAZ Brazil
VT+HF/ACCULUNG-UK United Kingdom


VT-Plus upgrades

(adds HF capability and RT-200 mode)

8831007 VT PLUS HF hardware and firmware factory service upgrade (for units lower than hardware v1.01.01; additional flat-rate charge required for factory service/calibration)

Standard accessories

9VT0015 Users Manual
8830200FG VT for Windows PC Software
75034 Serial Cable
1HD0011 Tilt Stand
 Power Cord (country specific)
VT-PLUS-7001 Accessory Kit (includes 16 accessories)



VT for Windows PC Software



VT PLUS HF standard accessories

Optional accessories

5022010 Soft Vinyl Carrying Case for VT PLUS HF

9530-0066 Hard-Sided Protective Carrying Case for VT PLUS HF (limited to stock on hand)

Test lungs

ACCU LUNG ACCU LUNG Portable Precision Test Lung (with Soft-Sided Carrying Case)

MI-14900 Michigan Instruments Non-Instrumented Single-Adult Test Lung

MI-11000 Michigan Instruments Non-Instrumented Dual-Adult Test Lung

MI-12952 Michigan Instruments Non-Instrumented Adult/Infant Test Lung

48499 Siemens 190 Test Lung

Parabolic airway resistors (for use with Michigan Instruments test lungs)

48129 Parabolic Airway Resistor ring

Printers

PRINTR/CTZ-US120V Printer 110 V, Citizen IDP 3110

PRINTR/CTZ-US220V Printer 220 V, Citizen IDP 3110

71072 Parallel Printer Cable, D25M-C36M

61096 Printer 120 V Power Supply

61097 Printer 220 V Power Supply

97116 DPU-414 and DPU-411 Printer Paper (minimum 7 rolls - price is per roll)

Accessory kit parts

1XX0015 Filter, External (Bacterial), 1 each

49343FG Adapter, DISS O2 Nut and Nipple with 1/4 in I.D. Hose Barb, 1 each

1FT0050 Tubing Adapter, Directional 15 mm OD x 15 mm OD), 2 each

1FT0049 Tubing Adapter (22 mm OD x 22 mm ID), 2 each

1FT0048 Tubing Adapter (22 mm OD x 22 mm OD), 2 each

1FT0045 Tubing Adapter (15 mm OD x 22 mm OD), 2 each

1FT0046 Tubing Adapter (15 mm OD x 15 mm OD), 2 each

1FT0047 Tubing Adapter (15 mm ID x 15 mm OD), 2 each

1FT0051 Tubing Adapter, Narrow Bore, 2 each

48478 Barb (Luer Lock - Male to 1/89 in ID tubing), 2 each

1FT0043 Tubing Adapter (1/4 in NPT Male to 1/8 in ID Tubing Barb Fitting), 2 each

1FT0005 Tubing Adapter (Luer Lock 1/16 in to Bulk-head Connection), 2 each

2FU0005 Fuse (500 mA)

67535 Tubing 1/8 in 4 ft long, 2 each

About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-0 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 and ISO 13485 medical device certified and our products are:

- CE Certified, where required
- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required
- NRC Compliant, where required

Fluke Biomedical.

Better products. More choices. One company.

Fluke Biomedical

6045 Cochran Road
Cleveland, OH 44139-3303 U.S.A.

Fluke Biomedical Europe

Science Park Eindhoven 5110
5692EC Son, The Netherlands

For more information, contact us:

In the U.S.A. (800) 850-4608 or
Fax (440) 349-2307
In Europe/M-East/Africa +31 40 267 5435 or
Fax +31 40 267 5436
From other countries +1 (440) 248-9300 or
Fax +1 (440) 349-2307
Email: sales@flukebiomedical.com
Web access: www.flukebiomedical.com

©2006-2011 Fluke Biomedical. Specifications subject to change without notice. Printed in U.S.A.
1/2011 2817925C D-EN-N

Modification of this document is not permitted without written permission from Fluke Corporation.