

DISCLAIMER FOR COMMUNICATION INTERFACES

Fluke provides communications interfaces (“as provided in the relevant communications interface document”) for development and use by its customers (“Users”) for their own use and within User’s controlled environment. Fluke uses and has verified the functionality of these communications interfaces in accordance with its own in-house performance standards using accepted test procedures.

Except for liability which cannot be excluded by law, Fluke shall have no responsibility for User’s development or use of these communications interfaces nor for any loss, corruption or limitation of access to the communications interfaces. Fluke is not responsible for any trouble shooting nor is it responsible for any damages incurred to any device used for testing (whether a Fluke device or third-party device), the equipment being tested or any person. Fluke will not be responsible for any injuries sustained due to unauthorized equipment modifications.

These communications interfaces are provided to User as-is and provided without warranty of any kind, whether statutory, written, oral, express or implied (including any warranties of merchantability or fitness for a particular purpose or any warranties arising from course of dealing or usage of trade). Fluke does not warrant that the communications interfaces will be delivered free of any interruptions, delays, omissions or errors or in a secure manner or that any faults or trouble shooting will be corrected.

Follow any other instructions in the relevant communications interface document and do not touch the Fluke device while deploying the communications interfaces as this may result in electrical shock hazards or improper operation.

No responsibility is assumed by Fluke for the use or reliability of communications interfaces that are not supplied by Fluke.

FLUKE Biomedical IDA-5 User Communication Interface

Revision 1.0 05/06/2020

The IDA-5 user virtual COM port over USB. Set Com port to 115200, No parity, 8 data, 1 Stop, no handshake,

GENERAL FORMAT

All commands have the following general format.

[COMMAND,param 1,param 2,...,param n]<CR><LF>

where the opening [and closing] are mandatory

IDA-5 User Communication Interface

COMMAND is one of the commands described in the table below.

,param 1 - ,param n are parameters required by some commands and must be separated by commas.

<CR> is the ASCII carriage return character.

<LF> is the ASCII line feed character.

If any command is not interpreted correctly by the IDA-5 the instrument will respond with [BADCMD]<CR><LF>.

COMMAND	Param number	Purpose	
BYE		Terminates computer control	NONE
POLL		Set the IDA into poling mode in which it will send data when requested by the host	A string that indicates the channels fitted as follows: [POLL,1,2,3,4] <CR><LF> Where each digit represents a channel. Functioning channels will be represented by the channel number. If a channel is not functioning then a 0 (zero) will replace the channel number: [POLL,1,2,0,4] <CR><LF> means that channels 1, 2, and 4 are functioning but channel 3 is not.
LOG		Set the IDA into logging mode in which it will sends data as it is generated by the measuring modules.	A string that indicates the channels fitted as follows: [LOG,1,2,3,4] <CR><LF> Where represents a functioning channel. If a channel is not functioning then a 0 (zero) will replace the channel number: [LOG,1,2,0,4] <CR><LF> means that channels 1, 2, and 4 are functioning but channel 3 is not. After this: as data becomes available, it is transmitted in the format described following this table.
CnF CnV C1FA		Start a flow test. where n is the channel number.	[OK] <CR><LF>
	1	Control Number.	
	2	Operator (name or initials).	
	3	Set Flow rate.	
CnO CnP		Start an occlusion test. where n is the channel number.	[OK] <CR><LF>
	1	Control Number.	
	2	Operator.	
	3	Set Flow Rate.	
CnPCA		Start a PCA test	[OK] <CR><LF>
	1	Control Number	
	2	Operator	
	3	Set Flow Rate	
	4 to 10	Reserved – not implemented	

IDA-5 User Communication Interface

COMMAND	Param number	Purpose	
END		End a test	[OK] <CR><LF>
	1	Channel number to end test on.	
FLOW		Request for current flow rate	[FLOW,nnnn.nn,hh:mm:ss.mmm]<CR><LF> where: nnn.nn Is the flow rate in ml/h. hh:mm:ss.mmm Time since start of test formats as: hh Hours mm Minutes ss Seconds mmm Milliseconds.
	1	The channel to retrieve the result from.	
VOL		Request for current Volume	[VOL,vvvv.vv,hh:mm:ss.mmm]<CR><LF> where: vvv.vv Is the volume in ml. hh:mm:ss.mmm Time since start of test formats as: hh Hours mm Minutes ss Seconds mmm Milliseconds
	1	The channel to retrieve the result from.	
PRES		Request for current pressure	[PRES,pppp,hh:mm:ss.mmm]<CR><LF> where: pppp Is the pressure in mmHg. hh:mm:ss.mmm Time since start of test formats as: hh Hours mm Minutes ss Seconds mmm Milliseconds
	1	The channel to retrieve the result from.	
RECS		Request number of records saved in flash memory.	[RECS,nn]<CR><LF> Where nn is the number of records (0 to 999) decimal.
GETREC		Request a record from the flash memory.	The header information followed by a series of strings in the same format as LOG mode.
	1	The record number to retrieve.	
DELREC		Delete single record from flash	
	1	The record number to delete.	
DELALL		Delete all records from flash.	[ERASED]<CR><LF> when finished.
STATUS		Request status string.	[STAT,ccccmmmmssss] where: c last command sent measuring module m current mode for channel s last status from measuring module

IDA-5 User Communication Interface

COMMAND	Param number	Purpose	
GETSN		Request serial number and firmware version	
	1	Channel selector 0 Main board 1 - 4 Measuring module	
GETHEAD		Request the report heading.	[HEAD,line 1,line2,line3]
SETHEAD		Send a new heading.	[OK] <CR><LF>
	1	Line 1 of the new heading.	
	2	Line 2 of the new heading.	
	3	Line 3 of the new heading.	
GETPARAMS		Get parameters and calibration factors.	
	1	Channel number to get parameters from. 0 Main instrument parameters 1 - 4 Measuring module parameters for the appropriate channel..	.
STATUS		To return the current status of all channels.	
GETTMPLTCOUNT		Get the number of templates saved in the IDA-5	
GETTMPLT		Retrieve a template from the IDA-5	
	1	The template number to retrieve	
SETTMPLT		Send template to IDA-5. This requires a sequence of SETTMPLT commands. e.g. [SETTMPLT,0,TemplateName,Comment] [SETTMPLT,n,Type,Rate,VolPress,unit,hh,mm,tol%] where n is template step number (1 to 6) type is FLOW, OCCL, PCA, BOLUS VOLPress, unit, hh, mm, tol% are values for step [SETTMPLT,END]	

FORMAT OF DATA IN LOG MODE.

When in logging mode data is sent to the host computer as it becomes available. The format of the data string is described below.

nftttttttt vvvvvvvv pppp <CR><LF>

Where

n is the channel number (zero based i.e. 0 to 3)

f is a status flag that can have the following values

- : (colon) normal result
- b bubble detected
- a air lock detected - the test must be restarted
- o over pressure (on occlusion test).

ttttttt elapsed time since test started in mili-seconds formatted as eight hexadecimal digits.

vvvvvvvv volume delivered since test started in 1000/ml formatted as eight hexadecimal digits.

pppp pressure in mmHg formatted as four hexadecimal digits. This is a two's complement signed value.

Additional characters are reserved fields.

<CR><LF> Terminator.

FORMAT OF RECORDS RETURNED BY GETREC,N

GETREC returns heading information followed by data records.

Heading records have the general format

Hn,value 1,value 2,..., value n<CR><LF>

Where:

H indicates this is a heading record

N is the heading number - followed by a comma

value 1 Heading values separated by commas.

|
value n

IDA-5 User Communication Interface

Heading Number	Parameter Number	Meaning
1	1	Serial number of IDA-5 followed by firmware version number.
2	1	Channel number test was performed on.
	2	Serial number of measuring module followed by firmware version number of measuring module.
3	1	Report heading line 1
4	1	Report heading line 2
5	1	Report heading line 3
6	1	Control Number
7	1	Operator
8	1	Manufacturer of device under test.
9	1	Model of device under test
10	1	Serial number of device under test
11	1	Location of device under test
12	1	Type of test (Flow, Occlusion, PCA)
13	1	Selected flow rate
14	1	Selected volume to deliver
15	1	Selected Basal Rate
16	1	Selected infusion volume (PCA)
17	1	Selected infusion time (PCA)
18	1	Start date of test
	2	Start time of test
19		reserved
20		reserved