

DISCLAIMER FOR COMMUNICATION INTERFACES

Fluke provides communications interfaces (“as provided in the relevant communications interface document”) for development and use by its customers (“Users”) for their own use and within User’s controlled environment. Fluke uses and has verified the functionality of these communications interfaces in accordance with its own in-house performance standards using accepted test procedures.

Except for liability which cannot be excluded by law, Fluke shall have no responsibility for User’s development or use of these communications interfaces nor for any loss, corruption or limitation of access to the communications interfaces. Fluke is not responsible for any trouble shooting nor is it responsible for any damages incurred to any device used for testing (whether a Fluke device or third-party device), the equipment being tested or any person. Fluke will not be responsible for any injuries sustained due to unauthorized equipment modifications.

These communications interfaces are provided to User as-is and provided without warranty of any kind, whether statutory, written, oral, express or implied (including any warranties of merchantability or fitness for a particular purpose or any warranties arising from course of dealing or usage of trade). Fluke does not warrant that the communications interfaces will be delivered free of any interruptions, delays, omissions or errors or in a secure manner or that any faults or trouble shooting will be corrected.

Follow any other instructions in the relevant communications interface document and do not touch the Fluke device while deploying the communications interfaces as this may result in electrical shock hazards or improper operation.

No responsibility is assumed by Fluke for the use or reliability of communications interfaces that are not supplied by Fluke.

FLUKE Biomedical INCUII User Communication Interface

Version 1.0 05/06/2020

COM PORT SETTINGS

Settings for the COM port should be made by the program that opens and uses the COM port such as a terminal emulation program (HyperTerminal, Tera Term or other). The settings in Device Manager are usually irrelevant because they are overridden by the controlling program.

The COM port should be set to:

- 115,200 baud
- No parity
- 8 data bits
- 1 stop bit

READING FORMATS

1. Readings do not have leading zeroes, except when preceding the decimal point.
2. All reading values are preceded by a letter that specifies the type of reading.
3. The table below summarizes the reading formats.

Test	Character preceding reading value	Maximum number of digits before decimal point	Number of decimal places	Example
Air Convection Temperature	Txx.xx,yy.yy,zz.zz,uu.uu, ww.ww	3	2	"T22.33,22.52,22.12,22.32,22.15" For channels that are not connected do not return value. For example, if ch2 and ch3 are not connected, the return would be : "T22.33,,,22.12,22.15"
Humidity	H	3	1	"H99.1"
Sound	S	3	2	"S45.3"
Airflow	A	1	2	"A1.41"
Conduction Temperature	Rxx.xx,yy.yy,zz.zz,uu.uu, ww.ww	2	2	"R22.33,22.52,22.12,22.32,22.15" For channels that are not connected do not return value. For example, if ch3 is not connected, the return would be : "R22.33,22.52,,22.12,22.15"
K-type temperature	K	2	2	K23.45
Skin temperature	N	2	2	N25.33

Reading Formats

COMMAND SPECIFICATIONS

Unless specified otherwise:

- Commands return *.

GENERAL COMMANDS

IDENT	Ask for the instrument identification and firmware version.
Legal modes:	All modes
Returns:	Model number: INCUII , followed by comma, followed by firmware version number, including build for example: " INCUII,1.00.06 "

SN	Ask for instrument serial number.
Legal modes:	All modes
Returns:	Product serial number (alpha-numeric up to 10 characters). If no serial number is defined, returns "none".

RESET	Reset the instrument.
Legal modes:	All modes
Returns:	Does not return *, once reset the normal power on response is sent.

LOCAL	Go to Local control mode.
Legal modes:	RMAIN
Returns:	LOCAL

REMOTE	Go to Remote control RMAIN mode.
Legal modes:	LOCAL
Returns:	RMAIN

QMODE	Query the mode.
Legal modes:	All modes
Returns:	The Remote mode mnemonic per table above.

MEASUREMENT COMMANDS

Legal in **RMAIN**, **CAL** and **DIAG** Remote modes, unless power up error condition is present and unhandled:

SMPRATE= <i>param</i>	Set sampling time.
<i>param</i>	sampling time: 20 sec to 120 seconds (increments of 10 seconds)
SNSGRP= <i>param1, param2, param3, param4, ..., all</i>	<p>Set measurement sensor group.</p> <p>T1-5: Air convection temperature sensors H: Humidity K: K-type thermocouple S: Sound N: Skin temperature sensor R1-5: radiant warmer sensor</p> <p>Example: SNSGRP=T1, T2, H, S By sending the command above, the unit identifies these tests as one group and based on the sampling rate of the tests it will send out data to the COM port based on this info. When START command is sent, the unit starts sending results as one result group where readings are separated by ,. For example, if sampling rate for the test is 20 seconds, after START is sent, the unit will send result packet every 20 seconds: 25.3,24.9,75.1,65.2 After 20 seconds 25.4, 24.6, 75.2,65.3 ... When ready to finish this test, send END</p>
<i>param</i>	T1, T2, T3, T4, T5, H, S, A, K, N, R1, R2, R3, R4, R5
START	Start readings from measurement sensor group and send to COM port.
END	End readings from measurement sensor group.
SETTUNIT= <i>param</i>	Set temperature unit.
<i>param</i>	Temperature unit: F or C
SETAFUNIT= <i>param</i>	Set air flow unit.
<i>param</i>	Air flow unit: FT or MT
QTUNIT	Query temperature unit.
Returns	Temperature unit: F or C
QAFUNIT	Query air flow unit.
Returns	FT or MT

TEMPERATURE COMMANDS

QATEMP = <i>ch1, ch2, ch3, ch4, ch5</i>	Select air convection temperature measurement.
<i>ch</i>	Temperature channel There are 5 different temperature sensor and the user should be able to retrieve data using this command. For example, QATEMP = 1 will return the temperature reading for channel 1. QATEMP = 2 will return the temperature reading for channel 2. And so on Also, QATEMP =1,2,3 will return the temperature readings for channel 1,2,3. In this case, 25.3,25.5, 25.2 will be returned. QATEMP =4,5 will return the temperature readings for channel 4 and 5. In this case, 25.6, 25.7 will be returned. And so on.
Returns	Air convection temperature measurement for each specified channel.

QCTEMP = <i>ch1, ch2, ch3, ch4, ch5</i>	Select conduction temperature measurement.
<i>ch</i>	Temperature channel
Returns	Conduction temperature reading for each specified channel.

HUMIDITY COMMANDS

QRHUM	Take humidity measurement.
Returns	Humidity reading.

SOUND COMMANDS

QSOUND	Take sound measurement.
Returns	Sound reading.

AIRFLOW COMMANDS

QAFLOW	Take air flow measurement.
Returns	Air flow reading.

SKIN TEMPERATURE

QSKTEMP	Take skin temperature measurement.
Returns	Skin temperature reading.