

Manual Supplement

Manual Title:	VT650/VT900A Users	Supplement Issue:	1
Part Number:	FBC - 0111	Date:	6/19
Print Date:	August 2018	Page Count:	1
Revision/Date:	1		

This supplement contains information necessary to ensure the accuracy of the above manual.

FLUKE®

Biomedical

Change #1, 234

On page 46, under Oxygen concentration, add note [1] to the Accuracy.

Accuracy^[1]..... ±1 % (VT900), ±2 % (VT650)

- [1] Allow minimum 45-minute warmup time for the Analyzer/oxygen sensor. After warmup, add 1 % to Oxygen accuracy if since the last calibration:
- More than 8 hours has elapsed
 - Ambient temperature changes more than ±2 °C
 - Ambient humidity changes more than ±5 % RH
 - Absolute sensor pressure changes more than ±1 %.

On page 48, table 15, replace row four and add new row:

Breath Rate	BPM	0.1 bpm	1 to 150 bpm	1 % ^[2]
Breath Rate	BPM	1 bpm	150 to 1500 bpm	2 % ^[3]

On page 49, table 15, add the following notes:

- [2] Breath reading using flow, pressure, and external trigger.
- [3] Breath reading using high frequency special test mode. External trigger cannot be used. Flow, pressure, and volume readings at high frequency are unspecified. Typically, these readings meet normal frequency specifications.