

Technical data

# ESA715 Electrical Safety Analyzer

## Our easiest-to-use electrical safety analyzer, ever.

Revolutionizing the field of electrical safety testing, the ESA715 Electrical Safety Analyzer offers an innovative approach to automated testing. Compact yet powerful, this groundbreaking device is an indispensable tool for healthcare technology management professionals. With its ability to perform from basic troubleshooting to comprehensive analysis, the ESA715 ensures uncompromised functionality in the field and in facilities. Adaptable to a range of global electrical safety standards, this all-in-one device seamlessly integrates a safety analyzer and patient simulator, redefining efficiency in medical equipment testing.



## Key features

- 5" Touchscreen with intuitive user interface
- Small, portable, ergonomic design with an integrated handle
- OneQA-enabled workflow automation eliminates need for laptop
- Wireless capability allows you to run test procedures in the field and automatically sync results
- Ten applied parts posts with innovative LED lighting guidance (ECG & Banana)
- Insulation resistance from 50V, 100V, 250V and 500V
- On-board automation with automated test sequences for compliance to key electrical safety standards NFPA 99 / AAMI ES1, AS/NZS 3551, IEC 62353 and parts of EN 50678 / EN 50699, and IEC 60601-1
- ECG waveform tests and dual-lead measurements combine the functionality of a simulator and safety analyzer in a single test tool
- Ability to upload pictures, instructions, and text to procedures
- Easy data entry through barcode scanner, external keyboard or on-screen keyboard
- 20 A at 120 V equipment current capability
- Replaceable mains fuses keep your unit in the field and out of the repair shop
- Rigorously tested for rugged field applications with CE and CSA in addition to Fluke-quality ruggedness for long-term reliability
- Enjoy precision, compliance, and total peace of mind through our comprehensive service and coverage program, eliminating unexpected costs with our Protect+ Program
- Global support network delivering prompt service and peace of mind to Fluke Biomedical customers worldwide



### Integrated OneQA workflow automation software allows you to:

- Improve productivity, efficiency and accuracy by automating the execution and analysis of tests
- Standardize testing procedures across all technicians – have everybody work like your senior technicians
- Ditch the laptop – execute OneQA right from the ESA715, keeping you mobile
- Reduce confusion around related to connections and testing with onboard guidance instructions
- Reduce onboarding time standardized test procedures and an intuitive user-friendly interface
- Build, organize and share reports and data on a centralized platform
- Achieve compliance with regulatory agencies with improved traceability
- Collaborate real-time with your team on procedures, data analysis and reporting
- Lower the cost of and improvement to Quality Assurance programs



[flukebiomedical.com/oneqa](https://flukebiomedical.com/oneqa)

## General Specifications

Safety standard compliance	IEC 61010-1: Overvoltage category II, pollution degree 2 IEC 61010-2-034: Measurement CAT II 300 V
Measurement equipment standard compliance	IEC 61557-16:2014, except IP40 for certain equipment outlets
Dimensions (w x d x h)	214 x 207 x 92 mm (8.4 x 8.1 x 3.6 inches)
Weight	1.3 kg (2.7 pounds)
Ingress protection	IP40 per IEC 60529, excluding equipment outlet
Operating temperature	0 to +35 °C (+50 to +95 °F)
Operating humidity	10-90 %, non-condensing
Storage temperature	-20 to +60 °C (-4 to +140 °F)
Storage humidity	5-95 %, non-condensing
Battery charging temperature	8 to +28 °C (+46 to +82 °F)
Altitude	100-127 V ac mains voltage and ≤5 150 V on input jacks: ≤ 5000 m 200-240 V ac mains voltage and ≤5 300 V on input jacks: ≤ 2000 m
Battery	Built-in rechargeable lithium-ion (<3600 mAh, 3.6 V)
Battery life	Up to 2 h
Connectivity	1 x USB-C for PC communication, 2 x USB-A for peripherals
Display	5-inch touchscreen
Data storage	> 10 000 measurements
Power (region dependent)	90-132 V ac, 20 A MAX, 47-63 Hz 90-132 V ac & 180-264 V ac, 10 A MAX, 47-63 Hz 90-132 V ac & 180-264 V ac, 15 A MAX, 47-63 Hz 90-132 V ac & 180-264 V ac, 16 A MAX, 47-63 Hz

## CSA approval

US version, 90-132 V ac, 20 A MAX, 47-63 Hz	CSA approved
NEMA 6-15 version, 90-132 V ac & 180-264 V ac, 15 A MAX, 47-63 Hz	CSA approved

## Electromagnetic Compatibility (EMC)

See users manual for more information.



## General Specifications (continued)

Measurement Specifications	
Mains voltage	
Range	0-264 V ac rms, limited by power specification
Accuracy	$\pm(2 \% + 0.2 \text{ V})$
Point to point voltage	
Range	0-300 V ac rms, limited by altitude specification
Accuracy	$\pm(2 \% + 0.2 \text{ V})$
Protective earth resistance & Point to point resistance	
Range	0-20 $\Omega$
Accuracy	$\pm(1 \% + 0.01 \Omega)$ at $\leq 2 \Omega$
	$\pm(1 \% + 0.1 \Omega)$ at $>2 \Omega$
Test current	Square wave of at least $\pm 200 \text{ mA}$ for $\leq 2 \Omega$
Open circuit voltage	Max $\pm 24 \text{ V}$
Equipment Current	
Range	0-20 A ac rms, limited by power specification
Accuracy	$\pm(5 \% + 0.05 \text{ A})$
Max duty cycle	0-10 A: continuous
	10-15 A: 7 min on/3 min off
	15-20 A: 5 min on/5 min off
<b>Note:</b> Equipment outlet will be turned off if max duty cycle is exceeded.	
Leakage Current	
Modes	ac + dc (true rms), ac only, dc only Mains on applied part leakage test: ac only
Patient load	AAMI ES1:1993 IEC 60601-1:2005 / IEC 62353:2014 See users manual for more information.
Crest factor	$\leq 2$
Range	0 $\mu\text{A}$ - 20 mA
Accuracy	dc & 20 Hz - 0.5 kHz: $\pm(1 \% + 1 \mu\text{A})$
	0.5 kHz - 50 kHz: $\pm(2.5 \% + 1 \mu\text{A})$
	50 kHz - 1 MHz: $\pm(5 \% + 1 \mu\text{A})$
Mains on applied part (Applies to: mains on applied part leakage, direct applied part leakage, alternative applied part leakage, alternative equipment leakage)	
Test voltage	Mains voltage $\pm 5 \%$
Current limit	1 mA $\pm 25 \%$ at 115 V for AAMI ES1, NFPA 99 3.5 mA $\pm 25 \%$ at 230 V for IEC 60601-1, IEC 62353 and EN 50678 / EN 50699 7.5 mA $\pm 25 \%$ at 230 V for AS/NZS 3551
Additional uncertainty	Up to $\pm 2 \mu\text{A}$ for 120 V, up to $\pm 4 \mu\text{A}$ for 230 V

**Note:** For alternative equipment leakage, alternative applied part leakage, direct equipment leakage, and direct applied part leakage tests, the leakage is scaled to nominal mains as per IEC 62353. Specified accuracy is not valid if (unscaled) current limit is exceeded.

## Measurement Specifications (continued)

Insulation Resistance		
Test voltage 250 V dc and 500 V dc	Range	0.1-100 MΩ
	Accuracy	±(2 % + 0.2 MΩ) at ≤ 10 MΩ ±(7.5 % + 0.2 MΩ) at >10 MΩ
	Test voltage accuracy	+20 % / -0 %
Test voltage 50 V dc and 100 V dc	Range	0.1-20 MΩ
	Accuracy	±(10 % + 0.2 MΩ)
	Test voltage accuracy	+30 % / -0 %
Short-circuit current	2 mA ± 0.25 mA	
Maximum load capacitance	2 μF	

ECG simulation		
Frequency accuracy	±2 %	
Amplitude accuracy	±5 %, for 2 Hz square wave	
Waveforms	ECG complex	30, 60, 120, 180 and 240 BPM
	Square wave	0.125 and 2 Hz, 50 % duty cycle
	Triangle wave	2 Hz
	Pulse, 63 ms	30 and 60 BPM
	Ventricular fibrillation	

Respiration simulation	
Rate	Apnea (0 BrPM) and 10-100 BrPM in 10 BrPM steps
Waveforms	Normal
Inspiration : expiration ratio	1:1
Impedance baseline	1000 Ω ± 5 % between leads
Impedance variation	1 ± 0.15 Ω
Respiration lead	LL or LA, user selectable

## Ordering information

### Models/descriptions

5577854	ESA715 Electrical Safety Analyzer-US, 120 V
5577724	ESA715 Electrical Safety Analyzer-01-France, 230 V
5577736	ESA715 Electrical Safety Analyzer-02-Europe, 230 V
5577749	ESA715 Electrical Safety Analyzer-03-Israel, 230 V
5577846	ESA715 Electrical Safety Analyzer-04-Denmark, 230 V
5577751	ESA715 Electrical Safety Analyzer-05-Australia, 230 V
5577760	ESA715 Electrical Safety Analyzer-06-UK, 230 V
5577772	ESA715 Electrical Safety Analyzer-07-Swiss, 230 V
5577785	ESA715 Electrical Safety Analyzer-08-Thailand, 230 V
5577797	ESA715 Electrical Safety Analyzer-09-Japan, 100 V
5577805	ESA715 Electrical Safety Analyzer-10-North America, 220 V
5577810	ESA715 Electrical Safety Analyzer-11-Brazil, 230 V
5577822	ESA715 Electrical Safety Analyzer-12-India, 230 V
5577831	ESA715 Electrical Safety Analyzer-13-South Africa, 230 V
6046090	ESA715 Electrical Safety Analyzer-14-NEMA615

### Standard accessories

5517611	Getting started guide
5517610	OneQA quick guide
5579600	USB C to USB A cable
2427138	Probe test banana jack 4MM tip black w/cap 175-290-001 EVZFTP74SW00 (USA/AUS/ISR/THAI/JAPAN)
2099044	Probe test banana jack 4MM tip red w/cap 175-290-003 EVZFTP74RT00 (USA/AUS/ISR/THAI/JAPAN)
4151242	USA/NEMA outlet to NBR14136 socket (BRAZIL)
3326842	Null post adapter
2248650	Carrying case
LINE CORD	Country-specific power cord

### Accessory kits (country specific)

3111008	USA/AUS/ISR Accessory Kit: test lead set, TP1 test probe set, AC285 alligator clip set (ESA T/L kit, USA)
3111024	EUR Accessory Kit: test lead set, TP74 test probe set, AC285 alligator clip set (ESA T/L kit, EUR)

### Optional accessories

6078258	Barcode scanner
6078184	USB Wi-Fi dongle



## About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-0 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.



## Fluke Biomedical regulatory commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 and ISO 13485 medical device certified and our products are:

- CE Certified, where required
- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required
- NRC Compliant, where required

### Fluke Biomedical

*We empower our everyday heroes to focus only on protecting lives.*

**Fluke Biomedical**  
6920 Seaway Boulevard  
Everett, WA 98203

**For more information, contact us at:**  
(800) 850-4608 or Fax (440) 349-2307  
sales@flukebiomedical.com  
flukebiomedical.com

©2025 Fluke Biomedical.  
Specifications subject to change without notice.  
5/2025 22884a-en

Modification of this document is not permitted without written permission from Fluke Corporation.