

Manual Supplement

Manual Title:	IDA-5 Users Manual	Supplement Issue:	3
Print Date:	April 2013	Issue Date:	9/19
Revision/Date:	2, 10/14	Page Count:	1

This supplement contains information necessary to ensure the accuracy of the above manual.
This manual is distributed as an electronic manual on the following CD-ROM:

CD Title:	IDA-5
CD Rev. & Date:	4/2013
CD PN:	3976006

Change #1, 214, 750

On page 33, in the **General Specifications** section, replace the **Weight, Electromagnetic Environment, Emissions Classification** and **FCC** with:
with:

Weight	IDA-5/4 3.17 kg (7.0 lb)
	IDA-5/2 2.81 kg (6.2 lb)
	IDA-5/1 2.58 kg (5.7 lb)
Electromagnetic Compatibility (EMC)	
International	IEC 61326-1: Basic Electromagnetic Environment; IEC 61326-2-1 CISPR 11: Group 1, Class A <i>Group 1: Equipment has intentionally generated and/or uses conductively-coupled radio frequency energy that is necessary for the internal function of the equipment itself.</i> <i>Class A: Equipment is suitable for use in all establishments other than domestic and those directly connected to a low-voltage power supply network that supplies buildings used for domestic purposes. There may be potential difficulties in ensuring electromagnetic compatibility in other environments due to conducted and radiated disturbances.</i> <i>Caution: This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.</i> <i>Emissions that exceed the levels required by CISPR 11 can occur when the equipment is connected to a test object.</i>
USA (FCC).....	47 CFR 15 subpart B. This product is considered an exempt device per clause 15.103.

Change #2, 698

On page 9, replace the **Printer** section with:

Use a printer that supports PCL-5 or PCL-5e printer-command language to print test results and reports from the Product.