

X-ray Test Equipment



© 2021.06 Unfors RaySafe G21217, Rev A

All rights are reserved. Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without the prior written consent of the copyright owner.

CONTENT

RaySafe X2..... 4

 RaySafe X2 Prestige System 4

 RaySafe X2 R/F & MAM System..... 7

 RaySafe X2 R/F System..... 8

 RaySafe X2 MAM System..... 8

 RaySafe X2 CT System 9

 RaySafe X2 Light System..... 9

 RaySafe X2 Survey System..... 10

 RaySafe X2 Volt System..... 10

 RaySafe X2 Sensors..... 11

 RaySafe X2 Upgrades..... 12

 RaySafe X2 Accessories & Spare Parts 12

RaySafe X2 Solo 17

 RaySafe X2 Solo R/F 17

 RaySafe X2 Solo Dent..... 19

 RaySafe X2 Solo Options & Upgrades 19

 RaySafe X2 Solo Accessories & Spare Parts..... 20

RaySafe ThinX..... 23

RaySafe DXR+ 24

Phantoms & Test Tools..... 25



 Phantoms 25



 CT 25



 R/F 26



 Test Tools..... 27





RaySafe X2





Part No / Model		RaySafe X2 Prestige System
5080633 / 1506013	RaySafe X2 Prestige R/F, M, C, L, mAs	
	Complete system in standard case: RaySafe X2 Base Unit mAs RaySafe X2 R/F Sensor RaySafe X2 MAM Sensor RaySafe X2 CT Sensor RaySafe X2 Light Sensor RaySafe X2 Flexi Stand RaySafe X2 Bluetooth Adapter RaySafe X2 Power Supply Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates	
5080640 / 1506014	RaySafe X2 Prestige R/F, M, C, L	
	Complete system in standard case: RaySafe X2 Base Unit RaySafe X2 R/F Sensor RaySafe X2 MAM Sensor RaySafe X2 CT Sensor RaySafe X2 Light Sensor RaySafe X2 Flexi Stand RaySafe X2 Bluetooth Adapter RaySafe X2 Power Supply Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates	



Part No / Model		RaySafe X2 Prestige System
5080745 / 1506035	RaySafe X2 Prestige R/F, M, C, S, mAs	
	Complete system in standard case: RaySafe X2 Base Unit mAs RaySafe X2 R/F Sensor RaySafe X2 MAM Sensor RaySafe X2 CT Sensor RaySafe X2 Survey Sensor RaySafe X2 Flexi Stand RaySafe X2 Bluetooth Adapter RaySafe X2 Power Supply Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates	
5080750 / 1506036	RaySafe X2 Prestige R/F, M, C, S	
	Complete system in standard case: RaySafe X2 Base Unit RaySafe X2 R/F Sensor RaySafe X2 MAM Sensor RaySafe X2 CT Sensor RaySafe X2 Survey Sensor RaySafe X2 Flexi Stand RaySafe X2 Bluetooth Adapter RaySafe X2 Power Supply Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates	



Part No / Model		RaySafe X2 Prestige System
5080761 / 1506037	RaySafe X2 Prestige R/F, M, L, S, mAs	
	Complete system in standard case: RaySafe X2 Base Unit mAs RaySafe X2 R/F Sensor RaySafe X2 MAM Sensor RaySafe X2 Light Sensor RaySafe X2 Survey Sensor RaySafe X2 Flexi Stand RaySafe X2 Bluetooth Adapter RaySafe X2 Power Supply Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates	
5080777 / 1506038	RaySafe X2 Prestige R/F, M, L, S	
	Complete system in standard case: RaySafe X2 Base Unit RaySafe X2 R/F Sensor RaySafe X2 MAM Sensor RaySafe X2 Light Sensor RaySafe X2 Survey Sensor RaySafe X2 Flexi Stand RaySafe X2 Bluetooth Adapter RaySafe X2 Power Supply Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates	







Part No / Model		RaySafe X2 R/F & MAM System
5080595 / 1506009	RaySafe X2 R/F, MAM, mAs	
	Complete system in standard case: RaySafe X2 Base Unit mAs RaySafe X2 R/F Sensor RaySafe X2 MAM Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates	
5080602 / 1506010	RaySafe X2 R/F, MAM	
	Complete system in standard case: RaySafe X2 Base Unit RaySafe X2 R/F Sensor RaySafe X2 MAM Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates	


Part No / Model		RaySafe X2 R/F System
5080517 / 1506001	RaySafe X2 R/F, mAs	
	<p>Complete system in standard case:</p> <p>RaySafe X2 Base Unit mAs</p> <p>RaySafe X2 R/F Sensor</p> <p>RaySafe X2 Flexi Stand</p> <p>RaySafe X2 Power Supply</p> <p>Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs</p> <p>Various documentation</p> <p>Calibration certificates</p>	
5080521 / 1506002	RaySafe X2 R/F	
	<p>Complete system in standard case:</p> <p>RaySafe X2 Base Unit</p> <p>RaySafe X2 R/F Sensor</p> <p>RaySafe X2 Flexi Stand</p> <p>RaySafe X2 Power Supply</p> <p>Cables: 2 m, 5 m, and 5 m ext USB</p> <p>Various documentation</p> <p>Calibration certificates</p>	
Part No / Model		RaySafe X2 MAM System
5080539 / 1506003	RaySafe X2 MAM, mAs	
	<p>Complete system in standard case:</p> <p>RaySafe X2 Base Unit mAs</p> <p>RaySafe X2 MAM Sensor</p> <p>RaySafe X2 Flexi Stand</p> <p>RaySafe X2 Power Supply</p> <p>Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs</p> <p>Various documentation</p> <p>Calibration certificates</p>	
5080542 / 1506004	RaySafe X2 MAM	
	<p>Complete system in standard case:</p> <p>RaySafe X2 Base Unit</p> <p>RaySafe X2 MAM Sensor</p> <p>RaySafe X2 Flexi Stand</p> <p>RaySafe X2 Power Supply</p> <p>Cables: 2 m, 5 m, and 5 m ext USB</p> <p>Various documentation</p> <p>Calibration certificates</p>	

Part No / Model		RaySafe X2 CT System
5080574 / 1506007	RaySafe X2 CT, mAs	
	Complete system in standard case: RaySafe X2 Base Unit mAs RaySafe X2 CT Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates	
5080588 / 1506008	RaySafe X2 CT	
	Complete system in standard case: RaySafe X2 Base Unit RaySafe X2 CT Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates	
Part No / Model		RaySafe X2 Light System
5080706 / 1506031	RaySafe X2 Light, mAs	
	Complete system in standard case: RaySafe X2 Base Unit mAs RaySafe X2 Light Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates	
5080714 / 1506032	RaySafe X2 Light	
	Complete system in standard case: RaySafe X2 Base Unit RaySafe X2 Light Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates	

Part No / Model RaySafe X2 Survey System		
5080723 / 1506033	RaySafe X2 Survey, mAs	
	Complete system in standard case: RaySafe X2 Base Unit mAs RaySafe X2 Survey Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates	
5080738 / 1506034	RaySafe X2 Survey	
	Complete system in standard case: RaySafe X2 Base Unit RaySafe X2 Survey Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates	

Part No / Model RaySafe X2 Volt System		
5307115 / 1506059	RaySafe X2 Volt mAs	
	Complete system in standard case: RaySafe X2 Base Unit mAs RaySafe X2 Volt Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates	
5307103 / 1506058	RaySafe X2 Volt	
	Complete system in standard case: RaySafe X2 Base Unit RaySafe X2 Volt Sensor RaySafe X2 Flexi Stand RaySafe X2 Power Supply Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates	






Part No / Model		RaySafe X2 Sensors
4559135 / 1252010	RaySafe X2 R/F Sensor	
	Multi-parameter sensor with stacked diodes (includes calibration certificates)	
4559147 / 1252020	RaySafe X2 MAM Sensor	
	Multi-parameter sensor with stacked diodes (includes calibration certificates). Dose/Dose Rate & HVL: Mo/Al, Mo/Cu, Mo/Mo, Mo/Rh, Rh/Ag, Rh/Al, Rh/Cu, Rh/Rh, W/Ag, W/Al, W/Cu, W/Rh, W/Ti. kVp: Mo/Cu, Mo/Mo, Mo/Rh, Rh/Ag, W/Ag, W/Al, W/Cu, W/Rh, W/Ti. Supports scanning mammography	
4559158 / 1252030	RaySafe X2 CT Sensor	
	100 mm pencil shaped CT sensor for dose, dose length product, dose rate and time. Includes calibration certificates	
4559164 / 1252040	RaySafe X2 Light Sensor	
	Luminance and illuminance sensor (includes calibration certificates)	
4673386 / 1252060	RaySafe X2 Survey Sensor	
	For performing leakage measurements and for measuring scattered radiation from X-ray systems. The sensitive survey sensor can also be used to measure low dose rates in the primary field. Includes calibration certificates	
5220340 / 1252070	RaySafe X2 Volt Sensor	
	The RaySafe X2 Volt sensor is used for measuring voltage from the output voltage test points on X-ray machines. Includes calibration certificates	

Part No / Model	RaySafe X2 Upgrades	
5092863 / 1916001	Upgrade to RaySafe X2 Base Unit with mAs	
	Upgrade of the RaySafe X2 Base Unit to enable invasive mA and mAs measurements	




Part No / Model	RaySafe X2 Accessories & Spare Parts	
4560024 / 1922052	RaySafe X2 Storm Case	
	Heavy duty, waterproof case with custom inlay for a complete RaySafe X2 system	
4560036 / 1922053	RaySafe X2 Storm Case System Option¹⁾	
	Upgrade from a standard case to a storm case	
4560013 / 1922050	RaySafe X2 Standard Alu Case	
	Lockable aluminium case with customized inlay	
4908594 / 1922093	RaySafe X2 Storm Case for GE Kit	
	Storm Case for a complete RaySafe X2 system with room for RaySafe X2 GE holders	
4908608 / 1922111	RaySafe X2 Storm Case for GE Kit System Option¹⁾	
	Upgrade to the RaySafe X2 Storm Case for GE kit when purchasing a RaySafe X2 system	

¹⁾ Only available when ordering a new system.

Part No / Model		RaySafe X2 Accessories & Spare Parts	
4559783 / 1902080	RaySafe X2 Flexi Stand		
	Flexible stand and holder for positioning RaySafe X2 sensors		
4559790 / 1902084	RaySafe X2 Vacuum Holder		
	Sensor holder with suction cup for various placements. Attached with a flexible arm and compatible with parts from the RaySafe X2 Flexi Stand		
4787978 / 1902093	RaySafe X2 Suction Cup Holder		
	Sensor holder for easy mounting and placement in various applications		
4559803 / 1902087	RaySafe X2 Scanning MAM Holder		
	Sensor holder for easy placement of the RaySafe X2 MAM Sensor when measuring on Scanning MAM, (e.g. Philips Microdose)		
4718087 / 1902091	RaySafe X2 Light Holder		
	Sensor holder for mounting the RaySafe X2 Light sensor on a monitor		



Part No / Model		RaySafe X2 Accessories & Spare Parts
4854002 / 1902088	RaySafe X2 Mammography Holder GE	
	Holder which fits the holder provided by GE. For the RaySafe X2 MAM Sensor on mammography machines from GE. Sensor position will be identical to ion-chamber setup	
4854016 / 1902089	RaySafe X2 Vascular Holder GE	
	Holder which adapts to the holder provided by GE. To be used with the RaySafe X2 R/F Sensor on angiography machines from GE. Sensor position will be identical to an ion-chamber setup	
4854025 / 1902090	RaySafe X2 Radiography Holder GE	
	Sensor holder which adapts to the holder provided by GE. To be used with the RaySafe X2 R/F Sensor on radiography machines from GE	
4853992 / 1922096	RaySafe X2 Power Supply Kit	
	International 12 V power supply including the most common power plugs	
4560097 / 3410034	USB Cable, 2 m, Green	
	USB to USB micro for sensor and PC connections	




Part No / Model RaySafe X2 Accessories & Spare Parts		
4560085 / 3410033	USB Cable, 5 m, Blue	
	USB to USB micro for sensor and PC connections	
4560106 / 3410035	USB Cable, 5 m Extender	
	USB male to USB female active extender to enable a 10 m long sensor cable	
4559674 / 1902026	mAs Cable, 2 m	
	Short mAs cable	
4559688 / 1902028	mAs Cable, 10 m	
	Long mAs cable	
4560060 / 1922064	RaySafe X2 Bluetooth Adapter	
	USB Bluetooth adapter for wireless communication to a PC	
4560072 / 1922068	RaySafe X2 GX Light Adapter	
	A GX adapter for Philips image intensifiers. To be used with the RaySafe X2 Light Sensor	

Part No / Model		RaySafe X2 Accessories & Spare Parts
4712697 / 1922067	RaySafe X2 MAM Filter, 2 mm Al	
	For kV measurements on Mo/Rh and W/Al above 40 kV	
5234544 / 1902103	RaySafe X2 Volt Cable - 10:1 Probe	
	The 10:1 probe attenuates the input signal by a factor of 10. Can be suitable for certain applications and machines e.g. Optima	
5234503 / 1902099	RaySafe X2 Volt Cable - Banana Connectors	
	For connection of the RaySafe X2 Volt Sensor to the voltage test points on X-ray machines	



RaySafe X2 Solo

RaySafe X2 Solo R/F

Part No / Model		RaySafe X2 Solo R/F Systems
5216000 / 1506053	RaySafe X2 Solo R/F with mAs, HVL, and Flexi Stand	
	RaySafe X2 Base Unit mAs RaySafe X2 R/F Sensor with active compensation: RaySafe X2 Power Supply RaySafe X2 Solo Case RaySafe X2 Flexi Stand Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates W/AI beam quality. Sensor has active compensation. Measures kVp, dose, rate, pulses, time, mA, mAs, and HVL	
5080853 / 1506047	RaySafe X2 Solo R/F with mAs and HVL	
	RaySafe X2 Base Unit mAs RaySafe X2 R/F Sensor with active compensation: RaySafe X2 Power Supply RaySafe X2 Solo Case Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates W/AI beam quality. Sensor has active compensation. Measures kVp, dose, rate, pulses, time, mA, mAs, and HVL	





Part No / Model		RaySafe X2 Solo R/F Systems
5264705 / 1506055	RaySafe X2 Solo R/F with HVL	
	RaySafe X2 Base Unit RaySafe X2 R/F Sensor with active compensation: RaySafe X2 Power Supply RaySafe X2 Solo Case Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates W/AI beam quality. Sensor has active compensation. Measures kVp, dose, rate, pulses, time, and HVL	
5080792 / 1506041	RaySafe X2 Solo R/F with mAs	
	RaySafe X2 Base Unit mAs RaySafe X2 R/F Sensor with active compensation: RaySafe X2 Power Supply RaySafe X2 Solo Case Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates W/AI beam quality. Sensor has active compensation. Measures: kVp, dose, rate, pulses, time, mA, and mAs	
5080809 / 1506042	RaySafe X2 Solo R/F	
	RaySafe X2 Base Unit RaySafe X2 R/F Sensor RaySafe X2 Power Supply RaySafe X2 Solo Case Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates W/AI beam quality. Sensor has active compensation. Measures: kVp, dose, rate, pulses, and time	

RaySafe X2 Solo Dent

Part No / Model	RaySafe X2 Solo Dent Systems	
5080811 / 1506043	RaySafe X2 Solo DENT with mAs	
	RaySafe X2 Base Unit mAs RaySafe X2 Dent Sensor RaySafe X2 Panoramic holder RaySafe X2 Power Supply RaySafe X2 Solo Case Cables: 2 m, 5 m, 5 m ext USB, and 10 m mAs Various documentation Calibration certificates W/AI beam quality. Measures kVp, dose, rate, pulses, time, mA, and mAs	
5080827 / 1506044	RaySafe X2 Solo DENT	
	RaySafe X2 Base Unit RaySafe X2 Dent Sensor RaySafe X2 Panoramic holder RaySafe X2 Power Supply RaySafe X2 Solo Case Cables: 2 m, 5 m, and 5 m ext USB Various documentation Calibration certificates W/AI beam quality. Measures kVp, dose, rate, pulses, and time	


RaySafe X2 Solo Options & Upgrades

Part No / Model	Options & Upgrades
5085445 / NA	RaySafe X2 Solo HVL Get HVL and Total Filtration simultaneously with the other parameters
4967755 / 1916019	RaySafe X2 Solo to RaySafe X2 Upgrade Upgrade of the RaySafe X2 Solo Base Unit to a full RaySafe X2 Base unit which can handle all types of RaySafe X2 sensors.

Part No / Model		RaySafe X2 Solo Accessories & Spare Parts	
4560024 / 1922052	RaySafe X2 Storm Case	Heavy duty, waterproof case with custom inlay for a complete RaySafe X2 system	
4560036 / 1922053	RaySafe X2 Storm Case System Option¹⁾	Upgrade from standard case to storm case	
4959576 / 1922115	RaySafe X2 Solo R/F / DENT Standard Alu Case	Lockable aluminium case with customized inlay with space for RaySafe X2 Panoramic Holder	
4560013 / 1922050	RaySafe X2 Standard Alu Case	Lockable aluminium case with customized inlay	
4908594 / 1922093	RaySafe X2 Storm Case for GE Kit	Storm Case for a complete RaySafe X2 system with room for RaySafe X2 GE holders	
4908608 / 1922111	RaySafe X2 Storm Case for GE Kit System Option¹⁾	Upgrade to the RaySafe X2 Storm Case for GE kit when purchasing a RaySafe X2 system	


¹⁾ Only available when ordering a new system.


Part No / Model		RaySafe X2 Solo Accessories & Spare Parts	
4559783 / 1902080	RaySafe X2 Flexi Stand	Flexible stand and holder for positioning RaySafe X2 sensors	
4959583 / 1902097	RaySafe X2 Panoramic Holder	Holder for accurate positioning of the RaySafe X2 DENT or R/F sensor on panoramic machines. Includes 10 gafchromic film slips (part 5600136)	
4560114 / 5600136	Gafchromic Film XR-QA2 - 10 pcs	To be used with the RaySafe X2 Panoramic holder	
4559790 / 1902084	RaySafe X2 Vacuum Holder	Sensor holder with suction cup for various placements. Attached with a flexible arm and compatible with parts from the RaySafe X2 Flexi Stand	
4787978 / 1902093	RaySafe X2 Suction Cup Holder	Sensor holder for easy mounting and placement in various applications	
4853992 / 1922096	RaySafe X2 Power Supply Kit	International 12 V power supply including the most common power plug	



Part No / Model		RaySafe X2 Solo Accessories & Spare Parts	
4560097 / 3410034	USB Cable, 2 m, Green	USB to USB micro for sensor and PC connections	
4560085 / 3410033	USB Cable, 5 m, Blue	USB to USB micro for sensor and PC connections	
4560106 / 3410035	USB Cable, 5 m Extender	USB male to USB female active extender to enable a 10 m long sensor cable	
4559674 / 1902026	mAs Cable, 2 m	Short mAs cable	
4559688 / 1902028	mAs Cable, 10 m	Long mAs cable	
4560060 / 1922064	RaySafe X2 Bluetooth Adapter	USB Bluetooth adapter for wireless communication to a PC	

RaySafe ThinX

Part No / Model	RaySafe ThinX	
4559199 / 1302024	RaySafe ThinX RAD	
	Fully automatic X-ray meter for radiographic applications. Measures kVp, dose, rate, time, pulses and HVL. Includes calibration certificate, user's manual and a soft protective case.	
4559208 / 1302025	RaySafe ThinX RAD kVp	
	Fully automatic X-ray meter for radiographic applications. Measures kVp, time, and pulses. Includes calibration certificate, user's manual and a soft protective case	
4559213 / 1302026	RaySafe ThinX RAD Dose	
	Fully automatic X-ray meter for radiographic applications. Measures dose, rate, time, pulses and HVL. Includes calibration certificate, user's manual and a soft protective case	
4559186 / 1302023	RaySafe ThinX Intra	
	Fully automatic X-ray meter for dental intraoral measurements. Measures kVp, dose, rate, time, pulses and HVL. Includes calibration certificate, user's manual and a soft protective case	
4588527 / 1302024-R	RaySafe ThinX RAD - R version	
	Fully automatic X-ray meter for radiographic applications. Measures kVp, dose, rate, time, pulses and HVL. Includes calibration certificate, user's manual and a soft protective case	
4588530 / 1302025-R	RaySafe ThinX RAD kVp - R version	
	Fully automatic X-ray meter for radiographic applications. Measures kVp, time, and pulses. Includes calibration certificate, user's manual and a soft protective case	


4559213 / 1302026-R	RaySafe ThinX RAD Dose - R	
	Fully automatic X-ray meter for radiographic applications. Measures dose, rate, time, pulses and HVL. Includes calibration certificate, user's manual and a soft protective case	
4588553 / 1302023-R	RaySafe ThinX Intra - R	
	Fully automatic X-ray meter for dental intraoral measurements. Measures kVp, dose, rate, time, pulses and HVL. Includes calibration certificate, user's manual and a soft protective case	




RaySafe DXR+

Part No / Model DXR+		
4558862 / 1151011	RaySafe DXR+ Direct X-ray ruler for X-ray/light field alignment and beam centering in R/F and mammography applications	







Phantoms & Test Tools



Phantoms

Part No / Model		CT
4717403 / 1921058	RaySafe Pro-CT Dose Phantom	
	Nested head, body and pediatric phantom with holes to place a CT sensor for CTDI measurements. Adult body is 32 cm, adult head/pediatric body is 16 cm, and pediatric head is 10 cm. Frosted PMMA with density 1.19 g/cm ³ . Delivered in a wheel carried Storm Case with inlay to fit a RaySafe X2 Base Unit, CT and R/F sensor	
4559815 / 1921046	Adult Head and Body CTDI Phantom	
	Body with diameter 32 cm and 4 holes for CTDI and head phantom with diameter 16 cm and 5 holes for CTDI. Frosted PMMA with density 1.19 g/cm ³ . Delivered in a wheel carried Storm Case with inlay to fit a RaySafe X2 Base Unit, CT and R/F sensor	
4712548 / 1921056	Head & Body CT Phantom, GE	
	Body with diameter 32 cm and 4 holes for CTDI and head phantom with diameter 16 cm and 5 holes for CTDI. Frosted PMMA with density 1.19 g/cm ³ . Includes a couch adapter for GE for easy mounting. Delivered in a wheel carried Storm Case with inlay to fit a RaySafe X2 Base Unit, CT and R/F sensor	

Part No / Model	R/F	
5102570 / 1921059	RaySafe Pro-Fluoro 150 Phantom <p>Made for acceptance and constancy tests of digital and analog radiography and fluoroscopy equipment. Suitable for routine quality checks on over-couch tubes, under-couch tubes and C-arms. Used for collimation/beam alignment, dynamic range, spatial and contrast resolution, homogeneity, and beam perpendicularity. Comes with holder for easy positioning on a table or chest wall. Complies with DIN 6868-150 and DIN 6868-4. Comes with a rugged case with space for a RaySafe X2 system</p>	
4717397 / 1922078	RaySafe Pro-Digi Radiography Phantom <p>Made for acceptance and constancy tests of digital radiography equipment. Used for collimation/beam alignment, dynamic range, spatial resolution, contrast resolution, homogeneity. Includes Al and Cu filter, alignment cone, bucky mounting device. Complies with DIN 6868-13, and DIN 6868-58. Delivered in a rugged case</p>	
4712527 / 1921052	RaySafe P 3D Fluoro Phantom <p>Is used to test the cone beam functionality of fluoroscopic X-ray systems. After carrying out image quality testing on 2D image acquisition with RaySafe Pro-Fluoro 150, this phantom is used to evaluate the 3D capability of 3D reconstruction module. Complies with DIN 6868-150. Delivered in a rugged case</p>	

Test Tools

Part No / Model	Test Tools	
4829300 / 1922103	RaySafe Pro-Slit Camera	
	Slit camera for measuring focal spot size on most machines. 0.01 mm wide and 10 mm long slit	
4829317 / 1922104	RaySafe Pro-Stand Full	
	Holder with adjustable height, 350-600 mm, to be used with the RaySafe Pro-Slit Camera or our pinholes. Can be used as a holder when doing manual HVL measurements. Can be tilted 10 degrees for accurate positioning on mammography machines. Includes a tool for easy alignment. Comes in a rugged case	
4829321 / 1922110	RaySafe Pro-Stand Basic	
	Holder with adjustable height, 350-600 mm, to be used with the RaySafe Pro-Slit Camera or our pinholes. Can be used as a holder when doing manual HVL measurements. Comes in a light-weight card board box	
3263272 / 07-501-2000	High-precision X-ray Test Pattern 07-501-2000	
	1.0-4.8 LP/mm, 0.1 mm Pb	
3263381 / 07-523-2000	High-precision X-ray Test Pattern 07-523-2000	
	0.5-5.0 LP/mm, 0.1 mm Pb	
3263453 / 07-535	High-precision X-ray Test Pattern 07-535	
	0.6-5.0 LP/mm, 0.05 mm Pb	

3263482 / 07-539	High-precision X-ray Test Pattern 07-539 1.5-20.0 LP/mm, 0.025 mm Pb	
3263803 / 07-633	X-ray Pinhole Assembly 07-633 0.010 mm (for focal spot size 0.5-1.0 mm)	
3263674 / 07-613	X-ray Pinhole Assembly 07-613 0.030 mm (for focal spot size < 1.0 mm)	
3263695 / 07-617	X-ray Pinhole Assembly 07-617 0.075 mm (for focal spot size 1.0-2.5 mm)	
3263663 / 07-611	X-ray Pinhole Assembly 07-611 0.100 mm (for focal spot size > 2.5 mm)	