

Impulse 4000 Plug-In

FLUKE.

Biomedical

The Ansur Impulse 4000 Plug-In is an add-on module for the Ansur Test Automation software application. The plug-in is used in conjunction with the Impulse 4000 Defibrillator Analyzer to provide test automation capabilities using the Ansur PC software.

Introduction

This read me file contains the latest information for Fluke Biomedical's Ansur Impulse 4000 Plug-In.

History

This chapter describes the changes made between each new released version.

Release 1.0.2

- Communication problem resolved related to USB to serial converter.
- BugFix in load settings for Pacer test.

Release 1.0.1

- Recompiled with the new version of a component compatible to windows Vista.
- Implemented Norwegian language settings to open custom setup tab and opening/saving templates.

Release 1.0.0

- V1.0.0 is the first release of the Impulse 4000 plug-in.

Language Changes

This chapter describes the changes made to the language file in each new version. The list includes only versions where the language file was changed.

Release 1.0.0

- Initial revision.

Fluke Biomedical
PO Box 9090
Everett, WA 98206-9090 USA

Fluke Biomedical France
30, rue Paul Claudel
91000 EVRY, France

For more information call:
U.S.A. (800) 648-7952 or Fax (425) 446-5629
Europe / M-East / Africa +33 1 60 78 88 99 or Fax +33 1 60 78 68 39
Other countries +1 (425) 446-6945 or Fax +1 (425) 446-5629
Internet: <http://www.flukebiomedical.com>

© 2008 Fluke Biomedical. All rights reserved.