






Manual Supplement

Manual Title:	601PRO XL Safety Sheet	Supplement Issue:	2
Part Number:	Web-Only	Issue Date:	11/16
Print Date:	July 2010	Page Count:	1
Revision/Date:			

This supplement contains information necessary to ensure the accuracy of the above manual.

Change #1

Add the following to the **Symbols** table:

	Conforms to relevant North American Safety Standards.
	Conforms to relevant South Korean EMC Standards.
	Conforms to relevant Australian Standards
	Conforms to European Union directives
CAT II	Measurement Category II is applicable to test and measuring circuits connected directly to utilization points (socket outlets and similar points) of the low-voltage MAINS installation.
	Consult user documentation.

Under the **Warnings** add:

- **Do not use the Product around explosive gas, vapor, or in damp or wet environments.**
- **Use the Product only as specified, or the protection supplied by the Product can be compromised.**

Change #2, 358

On page 2, add the following to the list of Warnings:

- **Connect an approved three-conductor mains power cord to a grounded power outlet.**