

Manual Supplement

Manual Title: 07-605-7777 Users
Print Date: August 2006
Revision/Date: 4


Supplement Issue: 1
Issue Date: 3/12
Page Count: 1

This supplement contains information necessary to ensure the accuracy of the above manual.

Change #1

On page 1-3, replace the entire page with:

EZ-CR-DIN Phantom w/Patches Model 07-605-7777



Low Contrast (Visible circles 1 = 5.6 %, 2 = 4 %, 3 = 2.8 %, 4 = 2.0 %, 5 = 1.6 %, 6 = 1.2 %)