

Instrument Security Procedures

Model: PS410

Product Name: PS410 patient simulator

Instrument Description:

- **The PS410 Patient Simulator is a compact, high-performance simulator for patient monitor testing**
- **This handheld device simulates a full range of cardiac rhythms and a wide variety of ECG conditions. It includes pacemaker simulation, 35 arrhythmia selections, and adult and pediatric normal-sinus rhythms.**

Memory Description: PS410 only outputs data and does not have memory for storage

Memory Cleaning Instructions: N/A