

Supplement

Title: ESA609 Safety Sht. Supplement Issue: 1
Print Date: February 2014 Issue Date: 11/14
Part Number: 4370092 Page Count: 1
Revision/Date: 1

This supplement contains information necessary
to ensure the accuracy of the above manual.

Change #1, 187

On page 5, under **Leakage Current**, replace the Accuracy with:

Accuracy

DC to 1 kHz $\pm(1\% \text{ of reading} + 1\ \mu\text{A})$

1 kHz to 100 kHz $\pm(2.5\% \text{ of reading} + 1\ \mu\text{A})$

100 kHz to 1 MHz $\pm(5\% \text{ of reading} + 1\ \mu\text{A})$