

# Manual Supplement

Manual Title: medTester 5000C Operators    Supplement Issue: **4**  
Part Number: 2243153    Issue Date: 4/18  
Print Date: February 2006    Page Count: 1  
Revision/Date:

---

---

This supplement contains information necessary to ensure the accuracy of the above manual.

## Change #1, 51108, 561, 571

On page 1-13, under **AAMI Load** replace the Accuracy section with the following:

- Accuracy  $\pm 1.0\%$  of reading DC + 3 LSD from 48 Hz to 1 kHz  
 $\pm 2.5\%$  of reading + 3 LSD from 1 kHz to 100 kHz  
 $\pm 5.0\%$  + 3 LSD from 100 kHz to 1 MHz