

Manual Supplement

Manual Title:	ProSim 4 Getting Started	Supplement Issue:	1
Part Number:	3931478	Date:	8/17
Print Date:	January 2011	Page Count:	1
Revision/Date:	3, 2/15		

This supplement contains information necessary to ensure the accuracy of the above manual.

Change #1, 217, 474

On page 4, replace the AC/DC Power Supply with:

AC/DC Power Supply	4766589
--------------------	---------

On page 16, under **General Specifications**, replace the Battery Charger

Battery Charger..... 110 Vac to 220 Vac, 50/60 Hz input, 6 V/3 A output. For best performance, the battery charger should be connected to a properly grounded ac receptacle.