

# SPOT Light

## SpO2 Functional Tester

### Addendum

### Introduction

This addendum gives the Instructions on how to switch between English and Chinese on units with Firmware version 1.03 or above. Use this addendum in conjunction with the SPOT Light Users manual.

The user may switch the language during power on by holding down the up-arrow button while pressing the SELECT button to power on the unit. This procedure shall change from the last language setting (English or simplified Chinese) to the new language (simplified Chinese or English). In order to switch back to the other language the unit must be powered off and restarted while holding the up arrow button.

### 附录

### 介绍

本附录作为对如何切换SPOT Light血氧功能测试仪的简体中文和英文语言界面的指导说明，适用于固件版本1.03或更高。本附录作为SPOT Light 使用说明书的补充与其共同使用。

用户可以通过以下操作切换中文和英文操作界面：在关机状态下，按住向上箭头键同时按下SELECT 按键至仪器开机。该操作将改变仪器上次开机运行时的语言界面（英文或简体中文）至另一个新语言界面（简体中文或英文）。如果需要切换回原来的语言，必须关机后重复上述操作。