



Manual Supplement

Manual Title:	VT Plus HF Operator	Supplement Issue:	1
Part Number:	2137275	Issue Date:	4/13
Print Date:	May 2005	Page Count:	1
Revision/Date:	2, 4/09		

This supplement contains information necessary to ensure the accuracy of the above manual. Enter the corrections in the manual if either one of the following conditions exist:

Change #1, 63303, 65379

On page 1-4, in the Symbols Table delete the **ETL** Symbol and add:

	Conforms to relevant North American Safety Standards.
	Conforms to relevant South Korean EMC Standards.

Under **Warnings** add:

- **Do not use the Product around explosive gas, vapor, or in damp or wet environments.**
- **Use the Product only as specified, or the protection supplied by the Product can be compromised.**

On page 1-11, add the following under **General Specifications:**

Safety IEC 61010-1: CAT II 300 V, Pollution Degree 2
 Operating Altitude..... 2000m
 Electromagnetic Environment IEC 61326-1: Basic
 Electromagnetic Compatibility Applies to use in Korea only. Class A Equipment (Industrial Broadcasting & Communication Equipment) ^[1]

[1] This product meets requirements for industrial (Class A) electromagnetic wave equipment and the seller or user should take notice of it. This equipment is intended for use in business environments and is not to be used in homes.