

Manual Supplement

Manual Title: SigmaPace 1000 Operators Supplement Issue: 3
Part Number: 2243306 Issue Date: 10/18
Print Date: July 2007 Page Count: 1
Revision/Date:

This supplement contains information necessary to ensure the accuracy of the above manual.

Change #1

On page 1-4, delete the following from the **Symbols** Table:

CE	Conforms to European Union directives
----	---------------------------------------

Change #2

On page 1-28, Table 1-3:

Change: 2184298

To: 4219453

On page 1-29, Table 1-6:

Change: 2184298

To: 4219453

Change #3, 632

On page 1-7, replace the third bullet with:

- A Sleep mode conserves the charge of the internal battery when not in use. Sleep mode is disabled when the external power supply is plugged into a suitable voltage source.