



# **Manual Supplement Rev. 1 for ULT800 TEE Transducer Leakage Current Tester (Rev. 2)**

Manual Title: ULT800 Users Guide	Supplement Issue: 1
Part Number: 2461434	Issue Date: 11/13
Print Date: August 2005	Page Count: 3
Revision/Date: May 2021	

---

---

This supplement contains information necessary to ensure the accuracy of the above manual.

## About This Manual Supplement

This document is a manual supplement against the ULT800 TEE Transducer Leakage Current Tester (Rev. 2 11/13). This document acts as an erratum to the original manual. Any changes detailed in this document reflect the most up-to-date information about the product.

## Change #1 CK 4/12/21

On page 6, under **Specifications** change:

From:

### Temperature

Operation.....15°C to 40°C

Storage.....15°C to 65°C

To:

### Temperature

Operation.....10°C to 40°C

Storage.....-20°C to 60°C