

Manual Supplement


Manual Title: VT305 Users
Print Date: January 2013
Revision/Date: 1

Supplement Issue: 3
Issue Date: 9/19
Page Count: 2

This supplement contains information necessary to ensure the accuracy of the above manual.

Change #1

On page 9, replace Figure 6 with:



gyo011.eps

Figure 6. Electrical Interfaces

Change #2, 610

On page 46, Table 7 add:

USB Cable	4034393
-----------	---------

Change #3, 238

On page 47, under **Dimensions** add:

Flow Port22 mm OD - ISO 5356-1:2015(E)

On page 47, under **Ambient pressure** add:

Altitude3000 m

On page 49, under **Oxygen concentration** add:

Sensor Technology.....Galvanic Fuel Cell

On page 50, under **Breathing Volume (V_{ti}, V_{te})** replace **Accuracy**, with:

Accuracy±2 %* or ±20 ml (> 6 slpm flow required)